

Product datasheet for **RC225217A1V**

Human CSRP3 (NM_001127656) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CSRP3 (NM_001127656) AAV Particle
Tag: Myc-DDK
Symbol: CSRP3
Synonyms: CLP; CMD1M; CMH12; CRP3; LMO4; MLP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225217 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCAACTGGGGCGGAGGCGCAAAATGTGGAGCCTGTGAAAAGACCGTCTACCATGCAGAAGAAATCC
AGTGCAATGGAAGGAGTTTCCACAAGACGTGTTTCCACTGCATGGCCTGCAGGAAGGCTCTTGACAGCAC
GACAGTCGCGGCTCATGAGTCGGAGATCTACTGCAAGGTGTGCTATGGGCGCAGATATGGCCCCAAAGGG
ATCGGGTATGGACAAGGCGCTGGCTGTCTCAGCACAGACACGGGCGAGCATCTCGGCCTGCAGTTCCAAC
AGTCCCCAAAGCCGGCACGCTCAGTTACCACCAGCAACCCTTCCAAATTCAGTGCAGTTTGGAGAGTC
CGAGAAGTGCCCTCGATGTGGCAAGTCAGTCTATGCTGCTGAGAAGGTTATGGGAGGTGGCAAGCCTTGG
CACAAGACCTGTTTCCGCTGTGCCATCTGTGGGAAGAGTCTGGAGTCCACAAATGTCACTGACAAAGATG
GGGAACCTTTATTGCAAGTTTGCTATGCCAAAAATTTGGCCCCACGGGTATTGGGTTTGGAGGCCTTAC
ACAACAAGTGGAAGAAGAAAGAA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225217 protein sequence
Red=Cloning site Green=Tags(s)

MPNWGGGAKCGACEKTVYHAEIQCNGRSFHKTCFHCMAKRALDSTTVAHSEIYCKVCYGRRYGPKG
IGYGQAGCLSTDTGEHLGLQFQQSPKPARSVTTSNPSKFTAKFGESEKPCRCGKSVYAAEKVMGGGKPW
HKTCFRCAICGKSLESTNVTDKGELYCKVCYAKNFGPTGIGFGLTQQVEKKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001127656
ORF Size:	582 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001127656.1, NP_001121128.1</u>
RefSeq Size:	1464 bp
RefSeq ORF:	584 bp
Locus ID:	8048
Cytogenetics:	11p15.1
MW:	21 kDa