

Product datasheet for **RC225232L3V**

CEP27 (HAUS2) (NM_001130447) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CEP27 (HAUS2) (NM_001130447) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CEP27
Synonyms:	C15orf25; CEP27; HsT17025
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130447
ORF Size:	612 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225232).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130447.1
RefSeq ORF:	615 bp
Locus ID:	55142
UniProt ID:	Q9NVX0
Cytogenetics:	15q15.2
MW:	23.1 kDa
Gene Summary:	The protein encoded by this gene is a subunit of the augmin complex. The augmin complex plays a role in microtubule attachment to the kinetochore and central spindle formation. [provided by RefSeq, Apr 2016]



[View online »](#)