

Product datasheet for **RC225259**

MFAP3 (NM_001135037) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MFAP3 (NM_001135037) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MFAP3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC225259 representing NM_001135037
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGTGTCTATTACATGATTGTTGCCTGATTGCCTTTACAATCACACTCATCTTGAATGTCACACGGC
 TGTGCATGATGAGCAGCCATCTTCGCAAGACTGAGAAGGCTATCAATGAGTTCCTTAGAACTGAAGGGGC
 TGAAACTTCAGAAGGCCTTTGAGATTGCAAAACGTATCCCCATCATTACCTCAGCCAAAACCTGGAG
 CTCGCCAAAGTCACACAATTTAAGACCATGGAGTTTGTCTGTTATATTGAAGAACTGGCAAGAAGTGCC
 CTCTCCACCTCTTATTCTAAACTGTCGAGCCTTTGTTGAGGAGATGTTGAGGCTGTGCGAGTGGACGA
 CCCTGATGACCTGGGTGAAAGAATTAAGAGAGACCTGCCTTGAATGCTCAAGGTGGCATCTATGTCATT
 AACCCAGAGATGGGACGGAGTAATTCACCAGGAGGAGATTCAGATGATGGCTCTCTGAATGAACAAGGCC
 AGGAAATAGCAGTTCAGGTTTCTGTCCACCTTCAGTCAGAAACCAAAAGTATTGATACAGAGTCTCAAGG
 CAGCAGTCATTTAGTCCACCTGATGATATAGGATCTGCAGAATCTAACTGTAACATAAAGATGGGGCA
 TATGAAAAGTGTGAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225259 representing NM_001135037
Red=Cloning site Green=Tags(s)

MSVYYMIVCLIAFTITLILNVTRLMMSSHLRKTEKAINEFFRTEGAEKLQKAFEIAKRIPITSAKLE
 LAKVTQFKTMEFARYIEELARSVPLPPLILNCRAFVEEMFEAVRVDDPDDLGERIKERPALNAQGGIYVI
 NPMEGRSNSPGGSDDGSLNEQGQEIYVQVSVHLQSETKSIDTESQSSHFSPDDIGSAESNCNYKDGA
 YENCQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001135037

ORF Size: 648 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001135037.2](#)

RefSeq Size: 4588 bp

RefSeq ORF: 651 bp

Locus ID: 4238
UniProt ID: [P55082](#)
Cytogenetics: 5q33.2
Protein Families: Transmembrane
MW: 24.1 kDa
Gene Summary: Component of the elastin-associated microfibrils.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for RC225259