

Product datasheet for **RC225261L4V**

SYCE2 (NM_001105578) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SYCE2 (NM_001105578) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SYCE2
Synonyms:	CESC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001105578
ORF Size:	654 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225261).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001105578.1
RefSeq ORF:	657 bp
Locus ID:	256126
UniProt ID:	Q6PIF2
Cytogenetics:	19p13.13
MW:	24.5 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is part of the synaptonemal complex formed between homologous chromosomes during meiotic prophase. The encoded protein associates with SYCP1 and SYCE1 and is found only where chromosome cores are synapsed. [provided by RefSeq, Dec 2012]