

Product datasheet for **RC225265A1V**

Human LINC02693 (NM_001113434) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LINC02693 (NM_001113434) AAV Particle
Tag: Myc-DDK
Symbol: LINC02693
Synonyms: FLJ12977; FLJ31874; FLJ33618; FLJ38937; FLJ41148
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225265 representing NM_001113434
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGAGAAAAATCACGACGGAAAGGCCAGCCCTCGCCATGCAGACGGGAAGCTTGGCCGCACGTGTG
ATCACCCCTTACGCCCTTGGAGTTTCACCCCCTCTTCGCGGGCTCCAACAGCCTGGGTACGGCTCCTTG
TCCTGTGTGGGCTCTCGGCTGCAAGAACACAGTCCCGAGCCCGGCGAGCCCGCGCCTCCCACCCGC
CGTGCACAGGCTGCACCTTACGCCCCAGCACTGCGCCTGCGGGACCCTGGATCGCTTTCCATTCTGA
TGACGTCTGTACTTCGTGGCTGCAGGCTCCTCAAGCCCCAGGTCTGTGCAGAGACGAGCAAAAGCAGCAG
AATATCCGTGCCACAGCTTTCTGGAGCACCCATTTTACTACCAGATTTGGAAGGAACGAAGTTAAGTAAT
TTCCAAGAATCATCTCCCTGCCCCACAACATGAAAGAAAAGACAAGAGAAGCACCCAGAAAGAGGAGG
GGAGAAGTGCCCCAGAAAAAATCATCCAGTCTCTGAAGCTTTGCCAGGTGGGCACAGGCCAGCCAGCCT
CTCCTCTGGTTGCCAGCTGGCTGCAGGCTGTCAATTAACCTTCGCCAGCATGCTTTTGAGTGTGCAG
AAGTGCTGCATGCCCTCTCTTTGAAAACCTGT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225265 representing NM_001113434
Red=Cloning site Green=Tags(s)

MGEKSRRKGPAPRHADGKLGRTCDHPYAPWSFTPSSRAPTAWVRPPCPVWASRLQEHSPPEPRRARAPPTR
RAQAALYAPALRLRDHLDRFSILMTSCTSWLQAPQAPGLCRDEQSSRISVPQLSGAPILLPDLEGTKLSN
FQESSPLPHKHERKDKRSTPEEEGRSAPEKIIQSLKLCPGGHRPASLSSGCPAGCRLSFNLPPSMMLSVQ
KCCMPSSLKTC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_001113434
ORF Size:	663 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001113434.2</u> , <u>NP_001106905.1</u>
RefSeq ORF:	666 bp
Locus ID:	339263
UniProt ID:	<u>A8MQB3</u>
Cytogenetics:	17p11.2
MW:	24.1 kDa