

## Product datasheet for RC225343L4V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## p53 (TP53) (NM\_001126115) Human Tagged ORF Clone Lentiviral Particle

**Product data:** 

**Product Type:** Lentiviral Particles

**Product Name:** p53 (TP53) (NM\_001126115) Human Tagged ORF Clone Lentiviral Particle

Symbol: TP53

Synonyms: BCC7; BMFS5; LFS1; P53; TRP53

**Mammalian Cell** 

Selection:

Puromycin

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

**ACCN:** NM\_001126115

ORF Size: 783 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC225343).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeq:** NM 001126115.1, NP 001119587.1

 RefSeq Size:
 2271 bp

 RefSeq ORF:
 786 bp

 Locus ID:
 7157

 UniProt ID:
 P04637

Cytogenetics: 17p13.1

**Protein Families:** Druggable Genome, Stem cell - Pluripotency, Transcription Factors





## p53 (TP53) (NM\_001126115) Human Tagged ORF Clone Lentiviral Particle - RC225343L4V

Protein Pathways: Amyotrophic lateral sclerosis (ALS), Apoptosis, Basal cell carcinoma, Bladder cancer, Cell cycle,

Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, Glioma, Huntington's disease, MAPK signaling pathway, Melanoma, Neurotrophin signaling pathway, Non-small cell lung cancer, p53 signaling pathway, Pancreatic cancer, Pathways in cancer, Prostate cancer,

Small cell lung cancer, Thyroid cancer, Wnt signaling pathway

MW: 29.6 kDa

Gene Summary: This gene encodes a tumor suppressor protein containing transcriptional activation, DNA

binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome.

Alternative splicing of this gene and the use of alternate promoters result in multiple

transcript variants and isoforms. Additional isoforms have also been shown to result from the

use of alternate translation initiation codons from identical transcript variants (PMIDs:

12032546, 20937277). [provided by RefSeq, Dec 2016]