

## Product datasheet for **RC225354L3V**

### **KCNK16 (NM\_001135107) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KCNK16 (NM_001135107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNK16
Synonyms:	K2p16.1; TALK-1; TALK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135107
ORF Size:	786 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001135107.1</a>
RefSeq ORF:	789 bp
Locus ID:	83795
UniProt ID:	<a href="#">Q96T55</a>
Cytogenetics:	6p21.2
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	28.6 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the family of potassium channel proteins containing two pore-forming P domains. This channel is an open rectifier which primarily passes outward current under physiological K<sup>+</sup> concentrations. This gene is expressed predominantly in the pancreas and is activated at alkaline pH. Several alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Sep 2008]