

Product datasheet for **RC225362A1V**

Human HCCS (NM_001122608) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HCCS (NM_001122608) AAV Particle
Tag:	Myc-DDK
Symbol:	HCCS
Synonyms:	CCHL; LSDMCA1; MCOPS7; MLS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225362 representing NM_001122608 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTTTGTCTCCATCTGCTCCTGCTGTTGCAGTTCAGGCCTCAAATGCTTCAGCGTCCCCACCTTCAG
GATGCCCGATGCATGAAGGAAAATGAAAGGCTGTCCAGTGAATACAGAGCCATCTGGCCCAACCTGTGA
GAAGAAACATACTCTGTGCTGCCACCGAGAACCGCCTATGAGTACGTGGAGTGTCCATTAGGGGC
ACTGCGGCTGAGAATAAGGAGAACCTAGATCCTTCAAATCTGATGCCACCACCAAATCAAACACCAGCTC
CAGATCAGCCATTTGCATTGTCTACTGTCAGAGAAGAGTCATCCATTCCGAGAGCAGATTAGAGAAAAA
GTGGGTTTACCCTTCTGAGCAGATGTTCTGGAATGCAATGTTAAAGAAAGGGTGAAGTGAAGGATGAG
GATATCAGTCAGAAGGATATGTATAATATCATTAGAATTCACAATCAGAATAACGAGCAGGCTTGAAGG
AGATTTTGAAGTGGGAAGCCCTTCATGCTGCAGAGTGCCTTGTGGTCCATCATTGATCCGGTTTGGAGG
GAAAGCAAAAGAGTATTCACCAAGGGCACGAATTCGTTCTGGATGGGGTATGAGTTGCCTTTTGATAGG
CACGATTGGATCATAAACCGTTGCGGGACAGAAGTTAGATATGTGATTGATTATTATGATGGTGGTGAAG
TCAACAAGGACTACCGATTCACCATCCTGGACGTCCTGCCTTAGATTCACTTTCGGCAGTATGGGA
CAGAATGAAAGTCGTTGGTGGCGTTGGACCTCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC225362 representing NM_001122608 Red=Cloning site Green=Tags(s) MGLSPSAPAVAVQASNASASPPSGCPMHEGKMKGCPVNTEPSGPTCEKKTYSVPAHQERAYEYVECPARG TAAENKENLDPSNLMPPPNQTPAPDQPFALSTVREESSIPRADSEKKWVYPSEQMFVNAMLKKGWKWKDE DISQKDMYNIIRIHNNQNEQAWKEILKWEALHAAECPCGPSLIRFGGKAKEYSPRARIRSWMGYELPFDR HDWIINRCGTEVRYVIDYYDGGEVNKDYQFTILDVRPALDSLAVWDRMKVAVWRWTS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001122608
ORF Size:	804 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001122608.1</u>
RefSeq ORF:	807 bp
Locus ID:	3052
UniProt ID:	<u>P53701</u>
Cytogenetics:	Xp22.2
MW:	30.4 kDa