

Product datasheet for **RC225376L3V**

TMEM72 (NM_001123376) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TMEM72
Synonyms:	C10orf127; KSP37
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001123376
ORF Size:	825 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225376).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001123376.1](#), [NP_001116848.1](#)

RefSeq ORF: 828 bp

Locus ID: 643236

UniProt ID: [A0PK05](#)

Cytogenetics: 10q11.21

Protein Families: Transmembrane



MW: 29.7 kDa

Gene Summary: This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep 2016]