

Product datasheet for **RC225378L4V**

CLEC12B (NM_001129998) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLEC12B (NM_001129998) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLEC12B
Synonyms:	UNQ5782
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001129998
ORF Size:	828 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225378).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001129998.1 , NP_001123470.1
RefSeq ORF:	831 bp
Locus ID:	387837
UniProt ID:	Q2HXU8
Cytogenetics:	12p13.2
Protein Families:	Transmembrane
MW:	31.4 kDa



[View online »](#)

Gene Summary:

Cell surface receptor that protects target cells against natural killer cell-mediated lysis. Modulates signaling cascades and mediates tyrosine phosphorylation of target MAP kinases. [UniProtKB/Swiss-Prot Function]