

## Product datasheet for **RC225382L4V**

### PPIL6 (NM\_00111298) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PPIL6 (NM_00111298) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPIL6
Synonyms:	bA425D10.6; dj919F19.1; PPIase; RSPH12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_00111298
ORF Size:	837 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225382).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_00111298.1</a> , <a href="#">NP_001104768.1</a>
RefSeq ORF:	1014 bp
Locus ID:	285755
UniProt ID:	<a href="#">Q8IXY8</a>
Cytogenetics:	6q21
MW:	31.3 kDa
Gene Summary:	Probable inactive PPIase with no peptidyl-prolyl cis-trans isomerase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)