

Product datasheet for **RC225388**

FRG2C (NM_001124759) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FRG2C
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC225388 representing NM_001124759
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGGGAAAGGGAATGAAGACCCCGATCTCCACTGTTCTCCATCCAGTGCTCCACTGACCAGCCCCCTT
TCCAACAGATCTCCTTTACAGAAAAGGGCTCAGATGAGAAGAAACCATTCAAAGGAAAAGGCAAGACCGC
CTTCTCCATTCCAGTGAGAAGCACACAAAGGCAAGCAGGATCAGATCCCAATCCAAACAAGGAGAAT
TCTGAGGAAACCAAGCTCAAGGCCGGGAACAGCACTGCTGGATCAGAACCAAGTCCAGCTCATATCAGG
AAACTGCAGGAAAAGAAAAATCAGTTCCAAGGATATCTGCCAAGACAGAGCAGGGAAGTGTCCAGAAGA
GGAGTGCAACTTGACGTTGAATAAAAAATCAAGATCCTCCACTGCTGTGCACAACAGTGAATCCAGGAG
ACCTGTGATGCCACCATAGGGGAAGTTCCAGGGCTTGCACTGGGCGCAGCAAGCGGCATAGGTCTCGGG
CCCTAGAAGTCCAAACACCGTCACTTCGAAAAAGCTTGGTGACCTCTGTGCGAGCTATGTCGGAGGCTGT
TTATCAAGACCTAGCCAGGTGTGGGCACAGCAGATCCATTCTCCACTGACCTGTGAGCAGCTGACACTG
CTCACTCGGCTCCGGGGCCTCTGTGTGCCAGGTGCAGACCTTGATTCCATGGCCACTCAGGCAGCTT
ATGTCTTCCTGCTGAGAGCTGGCTTGTCAGCCACACTGCCTGGTCTGGGGATTAGCCCTGGATAG
AGAAGCCCATCCCTTCCTGGGCAGGAGATAACTGAGCCTGTCACTGGATCAGATGAGGCTAAGCTGGGA
GCACCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC225388 representing NM_001124759

Red=Cloning site Green=Tags(s)

MGKGNEDPDLHCSSIQCSTDQPPFQQISFTEKGSDEKKPFKGKGTAFSHSSEKHTQRQAGSDPNPNKEN
SEETKLKAGNSTAGSEPESSSYQENCRKRKISSKDICQDRAGNCPEEECNLTNKKSRSSSTAVHNSEIQE
TCDAAHHRGSSRACTGRSKRHRRALEVQTPSLRKS LTVSVRAMSEAVYQDLAQVWAQQIHSPLTCEQLTL
LTRLRGPLCAQVQTLYSMATQAAYVFPAESWLVPATLPGPGDSALDREAHFPGQEITEPVSGSDEAKLG
AP

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8033_e06.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001124759

ORF Size:

846 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001124759.3

RefSeq ORF: 849 bp


Locus ID: 100288801

UniProt ID: A6NGY1


Cytogenetics: 3p12.3

MW: 30.6 kDa


Product images:



Circular map for RC225388



Western blot validation of overexpression lysate (Cat# [LY426628]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225388 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified FRG2C protein (Cat# [TP325388]). The protein was produced from HEK293T cells transfected with FRG2C cDNA clone (Cat# RC225388) using MegaTran 2.0 (Cat# [TT210002]).