

Product datasheet for **RC225391A1V**

Human KCTD15 (NM_001129994) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human KCTD15 (NM_001129994) AAV Particle
Tag: Myc-DDK
Symbol: KCTD15
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225391 representing NM_001129994
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCTCACCGCAAGGAGCGGCCGAGCGGGTCTCGCTTCACACACACGGCAGCACCGGCACCGCGGAGG
GAGGAAACATGTCCCGGCTGTCTCTACCCGGTCGCTGTGTCTCCCTGGCTGCCAGGGCATCCCCCT
GCCAGCCAGCTACCAAGTCCAATGCACCTGTGCACATCGATGTGGCGGCCACATGTACACCAGCAGC
CTGGCCACGCTACCAAGTACCCTGACTCCAGGATAAGCCGCCTCTTCAATGGCACTGAACCCATCGTCC
TGGACAGTTTGAAGCAACATTATTTATTGACCGGGATGGGAGATTTCCGCTACGTCCTGAGCTTCCT
GCGGACGTCCAAGCTGCTGCTTCCGGATGACTTAAAGGACTTCAGTCTGCTGTACGAGGAGGCGCGCTAC
TATCAGCTCCAGCCCATGGTGCAGAGCTGGAGCGCTGGCAGCAGGAGCAGGAGCAGCGCGCCGAGCC
GGGCCTGTGACTGCCTGGTGGTGCAGCTCACGCCGACTTGGCGAGCGGATCGCACTCAGCGCGAGAA
GGCCCTCATCGAGGAGGTCTTCCCCGAGACCGGAGACGTATGTGCAACTCCGTCAACGCCGGCTGGAAC
CAGGACCCACGCACGTATCCGCTTCCCGCTCAATGGCTACTGCCGGCTCAACTCGGTACAGGTCCTGG
AGCGGCTGTTCCAGAGGGGTTTCAGCGTGGCTGCGTCTGTGGGGCGGTGTGGACTCCTCCAGTTCAG
CGAGTATGTCTTTGCCGGGAGGAGCGGCCGCGCAGCCACCCCACTGCTGTTCAATCAAGCAGGAA
CCCCTGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC225391 representing NM_001129994 Red=Cloning site Green=Tags(s) MPHRKERPSGSSLHTHGSTGTAEGGNMSRLSLTRSPVSPLAAQGIPLPAQLTKSNAPVHIDVGGHMYTSS LATLTKYPDSRISRLFNGTEPIVLDSLKQHYFIDRDGEIFRYVLSFLRTSKLLLPDDFKDFSLLYEEARY YQLQPMVRELERWQQEQEQRRRSRACDCLVVRVTPDLGERIALSGEKALIEEVFPETGDVCMNSVNAGWN QDPHTVIRFPLNGYCRNLNSVQVLERLFQRGFSVAASCGGGVDSSQFSEYVLCREERRPQPTPTAVRIKQE PLD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001129994
ORF Size:	849 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001129994.1</u>
RefSeq ORF:	852 bp
Locus ID:	79047
UniProt ID:	<u>Q96S11</u>
Cytogenetics:	19q13.11
MW:	31.8 kDa