

Product datasheet for **RC225397L3V**

CTRP7 (C1QTNF7) (NM_001135171) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTRP7 (C1QTNF7) (NM_001135171) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1QTNF7
Synonyms:	CTRP7; ZACRP7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135171
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135171.1 , NP_001128643.1
RefSeq Size:	4477 bp
RefSeq ORF:	870 bp
Locus ID:	114905
UniProt ID:	Q9BXJ2
Cytogenetics:	4p15.32
Protein Families:	Secreted Protein
MW:	30.7 kDa



[View online »](#)