

## Product datasheet for RC225414L3V

### HYLS1 (NM\_001134793) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HYLS1 (NM_001134793) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HYLS1
Synonyms:	HLS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001134793
ORF Size:	858 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225414).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001134793.1</a>
RefSeq Size:	2073 bp
RefSeq ORF:	900 bp
Locus ID:	219844
UniProt ID:	<a href="#">Q96M11</a>
Cytogenetics:	11q24.2
MW:	32.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein localized to the cytoplasm. Mutations in this gene are associated with hydrolethalus syndrome. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Oct 2008]