

Product datasheet for **RC225429**

PEX26 (NM_001127649) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX26 (NM_001127649) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PEX26
Synonyms:	PBD7A; PBD7B; PEX26M1T; Pex26pM1T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC225429 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAAGAGCGATTCTTCGACCTCTGCAGCCCCCTCAGGGGGCTCGGGGACCCCTGCGCAGCAGCGAGC
CGGTGCGCGCGGTCCCGGCCGGGCGCCGGCGTGGACCTTCTGGAGGAGCGGCCGACCTCCTGGTGGT
GCACCTGGACTTCGGGCGGCGCTGGAGACCTGCGAGCGGGCTGGCAGAGTCTGGCCAACCACGCCGTG
GCAGAGGAACCCGCGGCACCTATTGGAGGTGAAGTGCTCCCTGTGTGTGTGGGGATCCAGGCCCTGG
CAGAAATGGATCGGTGGCAAGAAGTCTCTCCTGGGTCTTCAGTATTACCAGGTCCCTGAAAAGCTACC
CCCCAAAGTCTGGAGCTGTGCATTCTTTATACAGCAAAATGCAAGAGCCTGGAGCTGTGCTGGATGTG
GTGGGTGCCTGGCTCCAAGACCCAGCCAATCAAAACCTTCCAGAATATGGAGCCTTGGCAGAATTTACG
TGCAGCGGGTGCTGCTGCCTCTGGGCTGCTATCGGAGGCTGAGGAGCTAGTGGTGGGCTCTGCAGCCTT
TGGTGAGGAGCGGCGACTGGATGTACTTCAGGCCATTACACAGCGAGGCAGCAGCAGAAACAGGAACAC
TCAGGCTCTGAGGAGGCCAGAACCAACCTGGAAGGCTCTGTCTCCACAAGTTCTGTCACTACCGA
TGTTGGTTCGCCAGCTTTGGGACTCTGCGGTGAGCCACTTCTTTCTCTGCCCTTCAAAAAGAGTCTCCT
GGCTGCCTTGATCCTCTGTCTCCTGGTGGTGAGATTTGATCCAGCTTCCCTTCCCTCCCTGCACTTCCTC
TACAAGCTGGCCAGCTTCCGCTGGATCCGGAAGGCTGCATTTCTCGCCTCTACCAGCTCCGCATCC
GTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC225429 protein sequence
 Red=Cloning site Green=Tags(s)


MKSDSSTSAAPLRGLGGPLRSSEPVRAPAPAVDLLLEEAADLLVVHLDFAALETCEAWQSLANHAV
 AEEPAGTSLEVKCSL CVVGIQALAEEMDRWQEVLSWVLQYYQVPEKLPPKVLELCILLYSKMQEPGAVLDV
 VGAWLQDPANQNLPEYGALAEFHVQRVLLPLGCLSEAEELVVGSAAFGEERRLDVLQAIHTARQQQKEH
 SGSEEAQKPNLEGSVSHKFLSLPMLVRQLWDSAVSHFFSLPFFKKSLLAALILCLLVVRFDASPSSLHFL
 YKLAQLFRWIRKAAF SRLYQLRIRD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6307_h07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001127649

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001127649.3](#)

RefSeq Size: 4267 bp

RefSeq ORF: 918 bp

Locus ID: 55670


UniProt ID: [Q7Z412](#)

Cytogenetics: 22q11.21


MW: 33.9 kDa

Gene Summary: This gene belongs to the peroxin-26 gene family. It is probably required for protein import into peroxisomes. It anchors PEX1 and PEX6 to peroxisome membranes, possibly to form heteromeric AAA ATPase complexes required for the import of proteins into peroxisomes. Defects in this gene are the cause of peroxisome biogenesis disorder complementation group 8 (PBD-CG8). PBD refers to a group of peroxisomal disorders arising from a failure of protein import into the peroxisomal membrane or matrix. The PBD group is comprised of four disorders: Zellweger syndrome (ZWS), neonatal adrenoleukodystrophy (NALD), infantile Refsum disease (IRD), and classical rhizomelic chondrodysplasia punctata (RCDP). Alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Dec 2010]

Product images:



Circular map for RC225429



Western blot validation of overexpression lysate (Cat# [LY426835]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225429 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).