

Product datasheet for **RC225443A1V**

Human ASPA (NM_001128085) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ASPA (NM_001128085) AAV Particle
Tag:	Myc-DDK
Symbol:	ASPA
Synonyms:	ACY2; ASP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225443 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACTTCTTGTCACATTGCTGAAGAACATATACAAAAGGTTGCTATCTTTGGAGGAACCCATGGGAATG
AGCTAACCGGAGTATTTCTGGTTAAGCATTGGCTAGAGAATGGCGCTGAGATTCAGAGAACAGGGCTGGA
GGTAAAACCATTTATTACTAACCCAGAGCAGTGAAGAAGTGTACCAGATATATTGACTGTGACCTGAAT
CGCATTTTTGACCTTGAAAATCTTGGCAAAAAAATGTCAGAAGATTTGCCATATGAAGTGAGAAGGGCTC
AAGAAATAAATCATTTATTTGGTCCAAAAGACAGTGAAGATTCCTATGACATTATTTTTGACCTTCACAA
CACCACCTCTAACATGGGGTGCACCTTATTCTTGAGGATTCAGGAATAACTTTTTAATTCAGATGTTT
CATTACATTAAGACTTCTCTGGCTCCACTACCCTGCTACGTTTATCTGATTGAGCATCCTTCCCTCAAAT
ATGCGACCACTCGTTCCATAGCCAAGTATCCTGTGGGTATAGAAGTTGGTCCTCAGCCTCAAGGGTTCT
GAGAGCTGATATCTTGATCAAATGAGAAAAATGATTAACATGCTCTTGATTTTATACATCATTTCAAT
GAAGGAAAAGAATTTCTCCCTGCGCCATTGAGGTCTATAAAATTATAGAGAAAGTTGATTACCCCCGGG
ATGAAATGGAGAAATTGCTGCTATCATCCATCCTAATCTGCAGGATCAAGACTGGAAACCACTGCATCC
TGGGGATCCCATGTTTTTAATCTTGATGGGAAGACGATCCCACTGGGCGGAGACTGTACCGTGTACCCC
GTGTTTGTGAATGAGGCCGATATTACGAAAAGAAAGAAGCTTTTGCAAAGACAATAACTAACGCTCA
ATGCAAAAAGTATTCGCTGCTGTTTACAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC225443 protein sequence Red=Cloning site Green=Tags(s) MTSCHIAEEHIQKVAIFGGTHGNELTGVLVKHWLENGAEIQRTGLEVKPFITNPRAVKKCTRYIDCDLN RIFDLENLGKKMSEDLPYEVRRQAEINHLFGPKDSEDSYDIIFDLHNTTSNMGCTLIEDSRNNFLIQMF HYIKTSLAPLPCYVYLIEHPSLKYATTRSIKYPVGIEVGPQPQGVLRADILDQMRKMIKHALDFIHHFN EGKEFPPCAIEVYKIIIEKVDYPRDENGIEAAIIHPNLQDQDWKPLHPGDPMFLTLDGKTIPLGGDCTVYP VVFVNEAAYYEKKEAFAKTTKLTNNAKSIRCCLH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001128085
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001128085.1, NP_001121557.1</u>
RefSeq Size:	1368 bp
RefSeq ORF:	942 bp
Locus ID:	443
UniProt ID:	<u>P45381</u>
Cytogenetics:	17p13.2
MW:	35.7 kDa