

Product datasheet for **RC225452L3V**

Heme oxygenase 2 (HMOX2) (NM_001127205) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Heme oxygenase 2 (HMOX2) (NM_001127205) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Heme oxygenase 2
Synonyms:	HO-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001127205
ORF Size:	948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225452).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001127205.1 , NP_001120677.1
RefSeq Size:	1776 bp
RefSeq ORF:	951 bp
Locus ID:	3163
UniProt ID:	P30519
Cytogenetics:	16p13.3
Protein Families:	Transmembrane
Protein Pathways:	Porphyrin and chlorophyll metabolism



[View online »](#)

MW: 36 kDa

Gene Summary: Heme oxygenase, an essential enzyme in heme catabolism, cleaves heme to form biliverdin, which is subsequently converted to bilirubin by biliverdin reductase, and carbon monoxide, a putative neurotransmitter. Heme oxygenase activity is induced by its substrate heme and by various nonheme substances. Heme oxygenase occurs as 2 isozymes, an inducible heme oxygenase-1 and a constitutive heme oxygenase-2. HMOX1 and HMOX2 belong to the heme oxygenase family. Several alternatively spliced transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]