

Product datasheet for **RC225491A1V**

Human PCBP2 (NM_001128913) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PCBP2 (NM_001128913) AAV Particle
Tag:	Myc-DDK
Symbol:	PCBP2
Synonyms:	hnRNP-E2; HNRNPE2; HNRPE2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225491 representing NM_001128913 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGACACCGGTGTGATTGAAGGTGGATTAAATGTCACTCTCACCATCCGGCTACTTATGCATGGAAAGG
AAGTTGGCAGTATCATCGGAAAGAAAGGAGAATCAGTTAAGAAGATGCGCGAGGAGAGTGGTGACGAT
CAACATCTCAGAAGGGAATTGTCCTGAGAGAATTATCACTTTGGCTGGACCCACTAATGCCATCTTCAAA
GCCTTTGCTATGATCATTGACAACTGGAAGAGGACATAAGCAGCTCTATGACCAATAGCACAGCTGCCA
GTAGACCCCGGTCACCCTGAGGCTGGTGGTCCCTGCTAGTCAGTGTGGCTCTCTCATTGGAAAAGGTGG
ATGCAAGATCAAGGAAATACGAGAGAGTACAGGGGCTCAGGTCCAGGTGGCAGGGGATATGCTACCCAAC
TCAACTGAGCGGGCCATCACTATTGCTGGCATTCCACAATCCATCATTGAGTGTGTCAAACAGATCTGCG
TGGTCATGTTGGAGACTCTCTCCAGTCCCCCGAAGGGCGTGACCATCCCGTACCGGCCAAGCCGTC
CAGCTCTCCGGTCATCTTGCAGGTGGTCAGGCCTATACCATTCAAGGACAGTATGCCATTCCACAGCCA
GATTTGACCAAGCTGCACCAGTTGGCAATGCAACAGTCTCATTTCATGACGCATGGCAACACCGGAT
TCAGTGGCATTGAATCCAGCTCTCCAGAGGTGAAAGGCTATTGGGCAGGTTTGGATGCATCTGCTCAGAC
TACTTCTCATGAACTCACCATTCAAACGATTTGATTGGCTGCATAATCGGGCGTCAAGGCGCCAAAATC
AATGAGATCCGTCAGATGTCTGGGGCGCAGATCAAAATTGCGAACCCAGTGAAGGATCTACTGATAGGC
AGGTTACCATCACTGGATCTGCTGCCAGCATTAGCCTGGCTCAATATCTAATCAATGTCAGGCTTTCCTC
GGAGACGGGTGGCATGGGGAGCAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC225491 representing NM_001128913 Red=Cloning site Green=Tags(s) MDTGVIEGGLNVTLTIRLLMHGKEVGSIIIGKKGESVKKMREESGARINISEGNCPERIITLAGPTNAIFK AFAMIIDKLEEDISSSMTNSTAASRPPVTLRLVVPASQCGSLIGKGGCKIKEIRESTGAQVQVAGDMLPN STERAITTIAGIPQSIIECVKQICVVMLETLSQSPKGVITPYRKPSSSPVIFAGGQAYTIQGQYAIQP DLTKLHQLAMQQSHFPMTHGNTGFSGIESSSPEVKGYWAGLDASAQTTSHELTIPNDLIGCIIGRQGAKE NEIRQMSGAIKIANPVEGSTDRQVTITGSAASISLAQYLINVRLSSETGGMGSS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001128913
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001128913.1, NP_001122385.1</u>
RefSeq ORF:	1008 bp
Locus ID:	5094
UniProt ID:	<u>Q15366</u>
Cytogenetics:	12q13.13
MW:	35.2 kDa