

Product datasheet for **RC225528A1V**

Human CCRL2 (NM_001130910) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CCRL2 (NM_001130910) AAV Particle
Tag:	Myc-DDK
Symbol:	CCRL2
Synonyms:	ACKR5; CKRX; CRAM; CRAM-A; CRAM-B; HCR
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225528 representing NM_001130910 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCTACACCCGTTTCTTAAAAGGCAGTCTGAAGATGGCCAATTACACGCTGGCACCAGAGGATGAAT
ATGATGTCCTCATAGAAGGTGAAGTGGAGAGCGATGAGGCAGAGCAATGTGACAAGTATGACGCCAGGC
ACTCTCAGCCCAGCTGGTGCCATCACTCTGCTCTGTGTTTGTGATCGGTGTCCTGGACAATCTCCTG
GTTGTGCTTATCCTGGTAAAAATATAAAGGACTCAAACGCGTGAAAAATATCTATCTTCTAACTTGGCAG
TTTCTAACTTGTGTTTCTTGTACCTGCCCTTCTGGGCTCATGCTGGGGCGGATCCCATGTGTAAAAAT
TCTCATTGGACTGTACTTCGTGGGCTGTACAGTGAGACATTTTCAATTGCCTTCTGACTGTGCAAAAGG
TACCTAGTGTTTTTGCACAAGGGAACTTTTTCTCAGCCAGGAGGAGGGTGCCCTGTGGCATCATTACAA
GTGTCCTGGCATGGGTAAACAGCCATTCTGGCCACTTTGCCTGAATTCGTGGTTTATAAACCTCAGATGGA
AGACCAGAAATACAAGTGTGATTTAGCAGAACTCCCTTCCTGCCAGCTGATGAGACATTCTGGAAGCAT
TTTCTGACTTTAAAAATGAACATTTGCGTTCTGTCTCCCTATTTATTTTACATTTCTCTATGTGC
AAATGAGAAAAACACTAAGGTTCAAGGAGCAGAGGTATAGCCTTTTCAAGCTTGTTTTGCCATAATGGT
AGTCTTCTTCTGATGTGGGCGCCCTACAATATTGCATTTTCTGTCCACTTTCAAAGAACACTTCTCC
CTGAGTGACTGCAAGAGCAGCTACAATCTGGACAAAAGTTCACATCACTAACTCATGCCACCACCC
ACTGCTGCATCAACCCTCTCCTGTATGCGTTTCTTGATGGGACATTTAGCAAATACCTCTGCCGCTGTTT
CCATCTGCGTAGTAACACCCCACTTCAACCCAGGGGCGAGTCTGCACAAGGCACATCGAGGGAAGAACCT
GACCATTCACCGAAGTG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC225528 representing NM_001130910 Red=Cloning site Green=Tags(s) MIYTRFLKGSLKMANYTLAPEDEYDVLIEGELESDEAEQCDKYDAQALSAQLVPSLCSAVFVIGVLDNLL VVLILVKYKGLKRVENIYLLNLAVSNLCFLLTLPFWAHAGGDPCKILIGLYFVGLYSETFFNCLLTVQR YLVFLHKGNFFSARRRVPCGIITSVLAWVTAILATLPEFVVYKPMEDQKYKCAFSTPFLPADETFWKH FLTLMNISVLVPLFIPTFLYVQMRKTLRFREQRYSLFKLVFAIMVVFLMWAPYNIAFFLSTFKEHFS LSDCKSSYNLDKSVHITKLIATTHCCINPLLYAFLDGTFSKYLRCFHLRSNTPLQPRGQSAQGTSREEP DHSTEV SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001130910
ORF Size:	1068 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001130910.1</u>
RefSeq ORF:	1071 bp
Locus ID:	9034
UniProt ID:	<u>Q00421</u>
Cytogenetics:	3p21.31
MW:	40.8 kDa