

Product datasheet for RC225533

Dematin (DMTN) (NM_001114139) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dematin (DMTN) (NM_001114139) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dematin
Synonyms:	DMT; EPB49
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC225533 representing NM_001114139 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC

ATGGAACGGCTGCAGAAGGCCAAGATGGACAATCAGGTGCTGGGCTACAAGGACCTGGCTGCCATCCCCA
AGGACAAGGCCATCCTGGACATCGAGCGGCCGACCTCATGATCTACGAGCCTCACTTCACTTATCCCT
CCTGGAACACGTGGAGCTGCCTCGCAGCCGCGAGCGCTCGCTGTACCCAAATCCACATCCCCCCCACCA
TCCCCAGAGGTGTGGCGGACAGCCGGTCGCTGGAATCATCTCTCAGGCCTCGGCCCCAGAACCCTG
GAACCCCGGACCAGCCTGCCCATTTCCACCACCTGAGACCTCCCGCCAGATTCCAACATCTACAA
GAAGCCTCCCATCTATAAGCAGAGAGAGTCCGTGGGAGGCAGCCCTCAGACCAAGCACCTCATCGAGGAT
CTCATCATCGAGTCATCCAAGTTTCTGACGCCAGCCCGCCAGCCCAACCAGCCAGCCAAAATCGAAA
CCGACTACTGGCCATGCCCCCGTCTCTGGCTGTTGTGGAGACAGAATGGAGGAAGCGGAAGGCGTCTCG
GAGGGGAGCAGAGGAAGAGGAGGAGGAGGAAGATGACGACTCTGGAGAGGAGATGAAGGCTCTCAGGGAG
CGTCAGAGAGAGGAAGTCAAGGTTACTTCCAATTGGGAAAGATGATCTTGAAGAAGAGATGGAAA
AGTCATTGCCGATCCGAAGGAAAACCCGCTCTCTGCTGACCCGACACCCTTCCATACCTCCTTGCACCA
GGGAACGTCTAAATCTTCTCTCTCCCCGCTATGGCAGGACCACCCTGAGCCGGCTACAGTCCACAGAG
TTCAGCCCATCAGGGAGTGAGACTGGAAGCCCAGGCCTGCAGATCTATCCCTATGAAATGCTAGTGGTGA
CCAACAAGGGGCGAACCAAGCTGCCACCGGGGGTGGATCGGATCGGGCTTGAGAGGCATCTGTCTGCCGA
GGACTTCTCAAGGATTTTGCCATGTCCCTGAAGAGTTTGGCAAGCTGGCTCTGTGGAAGCGGAATGAG
CTCAAGAAGAAGGCCTCTCTCTTC

ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC225533 representing NM_001114139
Red=Cloning site Green=Tags(s)

MERLQKAKMDNQVLGYKDLAAIPKDKAILDIERPDLMIYEPHFYSLLEHVELPRSRERSLSPKSTSPPP
 SPEVWADSRSPGIISQASAPRTTGTPTSLPHFHPETSRLPDSNIYKPPPIYKQRESVGGSPQTKHLIED
 LIIESSKFPAAPDPNQPAKIEIDYWPCCPSLAVVETEWKRKASRRGAEEDDDSGEEMKALRE
 RQREELSKVTSNLGKMLKKEEMKSLPIRRKTRSLPDRTPFHTSLHQGTSKSSSLPAYGRITL SRLQSTE
 FSPSGSETGSPGLQIYPYEMLVVTNKGRTKLPVDRMLERHLSAEDFSRVFAMSPPEFGKLALWKRNE
 LKKKASLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8054_e08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN: NM_001114139

ORF Size: 1074 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001114139.4](#)

RefSeq ORF: 1077 bp

Locus ID: 2039


UniProt ID: [Q08495](#)

Cytogenetics: 8p21.3


MW: 40.5 kDa

Gene Summary: The protein encoded by this gene is an actin binding and bundling protein that plays a structural role in erythrocytes, by stabilizing and attaching the spectrin/actin cytoskeleton to the erythrocyte membrane in a phosphorylation-dependent manner. This protein contains a core domain in the N-terminus, and a headpiece domain in the C-terminus that binds F-actin. When purified from erythrocytes, this protein exists as a trimer composed of two 48 kDa polypeptides and a 52 kDa polypeptide. The different subunits arise from alternative splicing in the 3' coding region, where the headpiece domain is located. Disruption of this gene has been correlated with the autosomal dominant Marie Unna hereditary hypotrichosis disease, while loss of heterozygosity of this gene is thought to play a role in prostate cancer progression. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2014]

Product images:



Circular map for RC225533



Western blot validation of overexpression lysate (Cat# [LY426467]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225533 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).