

Product datasheet for **RC225549L3V**

SLC35G1 (NM_001134658) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SLC35G1
Synonyms:	C10orf60; POST; TMEM20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001134658
ORF Size:	1095 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225549).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001134658.1](#)

RefSeq ORF: 1098 bp

Locus ID: 159371

UniProt ID: [Q2M3R5](#)

Cytogenetics: 10q23.33

Protein Families: Transmembrane



MW: 39.7 kDa

Gene Summary: This gene encodes a transmembrane protein which is a member of the drug/metabolite transporter protein superfamily. The encoded protein may play a role in the regulation of calcium levels inside the cell. [provided by RefSeq, Sep 2016]