

Product datasheet for **RC225557L3V**

C6orf186 (METTL24) (NM_001123364) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C6orf186
Synonyms:	C6orf186
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001123364
ORF Size:	1098 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225557).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001123364.1](#), [NP_001116836.1](#)

RefSeq Size: 3152 bp

RefSeq ORF: 1101 bp

Locus ID: 728464

UniProt ID: [Q5JXM2](#)

Cytogenetics: 6q21



MW: 41.3 kDa