

Product datasheet for **RC225561L3V**

ZNF385B (NM_001113398) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF385B (NM_001113398) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF385B
Synonyms:	ZNF533
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113398
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113398.1 , NP_001106869.1
RefSeq Size:	2653 bp
RefSeq ORF:	1110 bp
Locus ID:	151126
UniProt ID:	Q569K4
Cytogenetics:	2q31.2-q31.3
MW:	50.4 kDa
Gene Summary:	May play a role in p53/TP53-mediated apoptosis.[UniProtKB/Swiss-Prot Function]


[View online »](#)