

Product datasheet for **RC225587**

Cyclin J (CCNJ) (NM_001134375) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Cyclin J
Synonyms:	bA690P14.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC225587 representing NM_001134375
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCTGGAGGGGCAGTGGTGGCGAGGACAGCTGGCCGCCGATATTCACCAAGCGCTTCGCTACAAGG
AGCTGAAGTTGCCCTCCTATAAAGGCCAGTCCCCTCAGTTAAGTCTCAGACGGTATTTTGCTGACTTGAT
TGCCATTGTGAGCAATCGCTTCACACTCTGCCCTTCTGCCGCCATCTTGCTGTCTATTTACTGGACCTG
TTTATGGACCGCTATGACATCTCTATCCAGCAGCTGCATTTAGTTGCGCTTTCCTGCCTGCTTCTAGCAA
GTAAATTTGAAGAAAAAGAACAGTGTGCCTAAGCTGGAGCAGCTCAACAGCCTGGGTTGTATGACTAA
TATGAATCTAGTATTAACAAAAACAAATTTGCTACATATGGAATATTATTATTAGAAACCTTTCAGTGG
AACCTCTGCCTTCAACAGCGCCCATTTTCATTGAGTATTATCTCTCTGAAGCAGTACACGAAACAGATC
TTCATGACGGCTGGCCAAATGATTTGCTTGGAAAAGACTAACTCTACATGGCCAAATATGCAGATTACTT
CCTGGAAGTATCTTTGCAAGATTATGCCTTTCTAAATTATGCACCTTCTTTAGTAGCTGCTGCATGTGTG
GCTTCTTCGAGGATTATACTTCGTCTTTCTCCAACGTGGCCTACAAGACTACATCGTCTTACTGCCTACT
CTTGGGATTTCTTAGTGCACTGATTGAACGACTGTTGATCGCTCATGATAATGATGTGAAAGAAGCAAA
CAAACAGAGAGGGCAAGCAGGACCTCAGTCAGCGCACTAAGTGTATTCAGACAGCCTCCAGCCATCA
CGGCCAGTTCACTTTAGCAACCTCAGTATCTCCATCAGACACATCAGACCTCACTGCAGTATCGGCATC
CTACGTCAGAACAACTGCTGTCAGCAGATTGTATCGACCACACACCTCATCTTACACACTACAGAC
ATGTCCTGCTGGCTTCCAACTAGTGTTCAGGGCCTTGGGCACATGCAGACTGGTGTGGGATGTCAGT
GCAATACAGTAGAAGTTAAGCCCTGTCTGAGTGTTCCTACAACCGGAGTTATCAGATAAATGAACATT
ACCCTTGATTAATCTCCATGTTTGAAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001134375.2](#)

RefSeq ORF: 1152 bp


Locus ID: 54619

UniProt ID: [Q5T5M9](#)


Cytogenetics: 10q24.1

MW: 43.7 kDa

Product images:



Circular map for RC225587



Western blot validation of overexpression lysate (Cat# [LY427405]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225587 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).