

Product datasheet for **RC225590A1V**

Human ZNF385 (ZNF385A) (NM_001130967) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	ZNF385
Synonyms:	HZF; RZF; ZFP385; ZNF385
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC225590 representing NM_001130967
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGATCCTCGGCAGCCTGAGCCGGGCAGGGCCCTGCCTCTGCTACGGCAGCCCCGATCATGCAGCCCC
CACTGGACCTCAAGCAGATCTGCCCTTCCCACTCGAGCCAGCCCCCTACCCTTGGCCTCTTCAGCACTA
CAGCACCATGGACCCTGTGCAGAAGGCTGTGCTCTCCACACTTTTGGGGGACCCTTGCTCAAGACCAAG
CGGCCCGTCATTTCTGTAATATCTGTCAAATCCGCTTCAATTCTCAGAGCCAGGCTGAGGCGCACTACA
AAGGTAATCGCCACGCCGACGAGTCAAAGGCATTGAGGCTGCCAAGACCAGAGGCAGGGAGCCTGGCGT
CCGAGAACCTGGAGACCCAGCTCCCCAGGCAGCACCCCAACAAATGGGGATGGTGTAGCACCCCGTCCA
GTTTCCATGGAGAATGGACTGGGGCCAGCCCCAGGATCCCCAGAGAAACAGCCTGGCTCCCCATCCCTC
CCAGCATTCCGGAGACTGGTCAGGGTGTAAACCAAGGGTGAAGGGGGGACTCCAGCCCCGGCTTCTTGCC
TGGGGGTAGCAAGGAAGAGGAGGAGAAAGCCAAGCGGCTGCTACTGTGCTCTGTGCAAGGTGGCTGTG
AACTCCCTGTCCAGCTTGAGGCACATAACAAAGGTAAGCACAAGACAATTCTGGAGGCCCGAAGTG
GGCTCGGGCCCATCAAAGCTTACCCTCGGCTGGGGCCTCCACCCCGGGGAACAGAGGCTCCTGCCCA
GGACCGAACTTTCACTGTGAGATCTGCAATGTCAAGGTCAACTCGGAGGTCCAAGTAAACAGCACATC
TCCAGCCGGCGGCACCGAGACGGCTGGCCGGGAAGCCCAACCCACTACTGAGCCGTCAAGAAGTCTA
GGGCGCCGGGGAGCTGGCGGCACGCTGACTTTCTCAAGGAGCTGCCAAGTCCCTGGCGGGCGGCT
GCTCCCCAGCCCCCTGGCGGTGGCTGCAGTGATGGCAGCGGCAGCAGGCTCGCCGCTGTCCCTGCCCGCG
GCTCCAGCCGACCTCTTCTCCAGGACCGCCGATCACTCACCTCTGCTTACCCGGCCCCCGGACCCA
TCCGAAGTGCACGACGACCATCCTCTTCTCCCGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATACGACGATAAGGTTTAA



Protein Sequence:	>RC225590 representing NM_001130967 Red=Cloning site Green=Tags(s) MILGSLSRAGPLLLRQPPIMQPPDLKQILPFPLEPAPTLGLFSNYSTMDPVQKAVLSHTFGGPLLKTK RPVISCNICQIRFNSQSQAIAHYKGNRHARRVKIEAAKTRGREPGVREPGDPAPPGSTPTNGDGVAPRP VSMENGLGPAPGSPEKQPGSPSPSIPETGQGVTKGEGGTPAPASLPGGSKEEEEKAKRLLYCALCKVAV NSLSQLEAHNKGTKHKITILEARSLGPIKAYPRLGPPTPGEPEAPAQDRTFHCEICNVKVNSEVQLKQHI SSRRHRDGVAGKPNLLSRHKKSRGAGELAGTLTFSELPKSLAGLLPSPLAVAAVMAAAAGSPLSLRP APAAPLLQGPPITHPLLHPAPGPIRTAHGPILFSY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001130967
ORF Size:	1158 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001130967.1</u>
RefSeq ORF:	1161 bp
Locus ID:	25946
UniProt ID:	<u>Q96PM9</u>
Cytogenetics:	12q13.13
MW:	40.3 kDa