

Product datasheet for **RC225590L3V**

ZNF385 (ZNF385A) (NM_001130967) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF385 (ZNF385A) (NM_001130967) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF385
Synonyms:	HZF; RZF; ZFP385; ZNF385
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130967
ORF Size:	1158 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225590).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130967.1
RefSeq ORF:	1161 bp
Locus ID:	25946
UniProt ID:	Q96PM9
Cytogenetics:	12q13.13
Protein Families:	Transcription Factors
MW:	40.3 kDa



[View online »](#)

Gene Summary:

Zinc finger proteins, such as ZNF385A, are regulatory proteins that act as transcription factors, bind single- or double-stranded RNA, or interact with other proteins (Sharma et al., 2004 [PubMed 15527981]).[supplied by OMIM, Oct 2008]