

Product datasheet for **RC225594L3V**

ZSCAN23 (NM_001012455) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZSCAN23
Synonyms:	dJ29K1.3; dJ29K1.3.1; ZNF390; ZNF453
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012455
ORF Size:	1167 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225594).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001012455.1
RefSeq Size:	3178 bp
RefSeq ORF:	1170 bp
Locus ID:	222696
UniProt ID:	Q3MJ62
Cytogenetics:	6p22.1



MW: 45 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]