

Product datasheet for RC225598

PTGER3 (NM_001126044) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTGER3 (NM_001126044) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PTGER3
Synonyms:	EP3; EP3-I; EP3-II; EP3-III; EP3-IV; EP3-VI; EP3e; Inc003875; PGE2-R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC225598 representing NM_001126044 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGAGACCCGGGGCTACGGAGGGGATGCCCCCTTCTGCACCCGCCTCAACCACTCCTACACAGGCA
TGTGGGCGCCGAGCGTTCCGCCGAGGCGGGGCAACCTCACGCGCCCTCCAGGGTCTGGCGAGGATTG
CGGATCGGTGTCCTGGCCTTCCCGATCACCATGCTGCTCACTGGTTTCGTGGCAACGCACTGGCCATG
CTGCTCGTGTGCGCAGCTACCGCGCCGGAGAGCAAGCACAAGAGTCTTCTGCTGTGCATCGGCT
GGCTGGCGCTCACCGACCTGGTCGGCAGCTTCTCACCACCCGGTCGTCATCGTCGTGTACCTGTCCAA
GCAGCGTTGGGAGCACATCGACCCGTCGGGGCGGCTCTGCACCTTTTCGGGCTGACCATGACTGTTTTTC
GGGCTCTCCTCGTTGTTTCATCGCCAGCGCCATGGCCGTCGAGCGGGCGCTGGCCATCAGGGCGCCGCACT
GGTATGCGAGCCACATGAAGACGCGTGCCACCCGCGTGTGCTGCTCGGCGTGTGGCTGGCCGTCGCTCGC
CTTCGCCCTGCTGCCGCTGCTGGGCGTGGGCCAGTACACCGTCCAGTGGCCCGGGACGTGGTGTTCATC
AGCACCGGGCGAGGGGGCAACGGGACTAGCTCTTCGCATAACTGGGGCAACCTTTTCTTCGCTCTGCT
TTGCCCTTCTGGGCTCTGGGCTGACAGTACCTTTTCTGCAACCTGGCCACCATTAAGGCCCTGGT
GTCCCGCTGCCGGGCAAGGCCACGGCATCTCAGTCCAGTGGCCAGTGGGCGCCATCACGACCGAGACG
GCCATTCAGCTTATGGGATCATGTGCGTGTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCT
AAATGATCTTCAATCAGACATCAGTTGAGCACTGCAAGACACACGAGAAAGCAGAAAGAAATGCAACTT
CTTCTTAATAGCTGTTCCGCTGGCTTCACTGAACCAGATCTGGATCCTTGGGTTTACCTGCTGTTAAGA
AAGATCTTCTCGAAAGTTTGGCAGATCAGGTACCACAAAACAATATGCATCCAGCTCCACCTCCT
TACCCTGCCAGTGTCTCAACCTTGATGTGGAGCGACCATTTGGAAAGA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC225598 representing NM_001126044
Red=Cloning site Green=Tags(s)


MKETRGYGGDAPFCTRLNHSYTGMPERSAEARGNLTRPPGSGEDCGSVSVAFPIITMLLTGFVGNALAM
 LLVSRSYRRRESKRKKSFLLCIGWLALTDLVGQLLTTPVVIVVYLSKQRWEHIDPSGRLCTFFGLTMTVF
 GLSSLFIASAMAVERALAIRAPHWYASHMKTRATRAVLGVWLVAFALLPVLGVGQYTVQWPWTWCFI
 STGRGGNGTSSSHNWGNLFFASAFALGLLALTVTFSCNLATIKALVSRCAKATASQSSAQWGRITTET
 AIQLMGIMCVLSVCWSPLLIMMLKMIFNQTSVEHCKTHTEKQECNFLLIIVRLASLNQILDWPVYLLLR
 KILLRKFQCQIRYHTNNYASSSTSLPCQCSSTLMWSDHLER

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8056_e10.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001126044

ORF Size: 1170 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001126044.2](#)

RefSeq Size: 2055 bp

RefSeq ORF: 1173 bp

Locus ID: 5733

UniProt ID: [P43115](#)

Cytogenetics: 1p31.1


Protein Families: Druggable Genome, GPCR, Transcription Factors, Transmembrane

Protein Pathways: Calcium signaling pathway, Neuroactive ligand-receptor interaction

MW: 43.3 kDa

Gene Summary: The protein encoded by this gene is a member of the G-protein coupled receptor family. This protein is one of four receptors identified for prostaglandin E2 (PGE2). This receptor may have many biological functions, which involve digestion, nervous system, kidney reabsorption, and uterine contraction activities. Studies of the mouse counterpart suggest that this receptor may also mediate adrenocorticotrophic hormone response as well as fever generation in response to exogenous and endogenous stimuli. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]

Product images:



Circular map for RC225598