

Product datasheet for **RC225608L4V**

PIP5KL1 (NM_001135219) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	PIP5KL1
Synonyms:	PIPKH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001135219
ORF Size:	1182 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225608).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001135219.1](#)

RefSeq ORF: 1185 bp

Locus ID: 138429

UniProt ID: [Q5T9C9](#)

Cytogenetics: 9q34.11

Protein Families: Druggable Genome



MW: 44.4 kDa

Gene Summary: PIP5KL1 is a phosphoinositide kinase-like protein that lacks intrinsic lipid kinase activity but associates with type I PIPKs (see PIP5K1A; MIM 603275) and may play a role in localization of PIPK activity (Chang et al., 2004 [PubMed 14701839]).[supplied by OMIM, Jun 2009]