

Product datasheet for **RC225637L4V**

METT5D1 (METTL15) (NM_001113528) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	METT5D1
Synonyms:	METT5D1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001113528
ORF Size:	1221 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225637).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001113528.1
RefSeq ORF:	1224 bp
Locus ID:	196074
UniProt ID:	A6NJ78
Cytogenetics:	11p14.1
MW:	45.9 kDa



Gene Summary:

Probable S-adenosyl-L-methionine-dependent methyltransferase.[UniProtKB/Swiss-Prot Function]