

Product datasheet for RC225640

Neuroserpin (SERPINI1) (NM_001122752) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Neuroserpin (SERPINI1) (NM_001122752) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SERPINI1
Synonyms:	HNS-S1; HNS-S2; neuroserpin; P112
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC225640 representing NM_001122752 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTTTCCTTGGACTCTTCTCTTTGCTGGTTCTGCAAAGTATGGCTACAGGGGCCACTTCCCTGAGG
AAGCCATTGCTGACTTGTCAAGTAAATGTATAATCGTCTTAGAGCCACTGGTGAAGATGAAAATATTCT
CTTCTCTCCATTGAGTATTGCTCTTGAATGGGAATGATGGAACCTGGGGCCCAAGGATCTACCCAGGAA
GAAATCCGCCACTCAATGGGATATGACAGCCTAAAAATGGTGAAGAATTTTCTTCTTGAAGGAGTTTT
CAAACATGGTAAGTCTAAAGAGAGCCAATATGTGATGAAAATTGCCAATTCCTTGTTGTGCAAAATGG
ATTTCAATGTAAGGAGTTTTTGAATGATGAAAAATATTTAATGCAGCAGTAAATCATGTGGAC
TTCAGTCAAAATGTAGCCGTGGCCAACACATCAATAAGTGGGTGGAGAATAACACAAACAATCTGGTGA
AAGATTTGGTATCCCAAGGGATTTTGTGCTGCCACTTATCTGGCCCTCATTAAATGCTGTCTATTTCAA
GGGGAAGTGAAGTGCAGTTTAGGCCTGAAAACTAGAACCTTTTCTTCACTAAAGATGATGAAAGT
GAAGTCCAAATCCAATGATGTATCAGCAAGGAGAATTTTATTATGGGAATTTAGTGTGGCTCCAATG
AAGCTGGTGGTATCTACCAAGTCTAGAAATACCATATGAAGGAGATGAAATAAGCATGATGCTGGTGTCT
GTCCAGACAGGAAGTTCCTCTTGTACTCTGGAGCCATTAGTCAAAGCACAGCTGGTTGAAGAATGGGCA
AACTCTGTGAAGAAGCAAAAAGTGAAGTATACCTGCCAGGTTACAGTGAACAGGAAATGATTTAA
AAGATGTTTTGAAGGCTCTTGAATAACTGAAATTTTCATCAAAGATGCAAAATTTGACAGGCCTCTATGA
TAATAAGGAGATTTTCTTTCCAAAGCAATTCACAAGTCTTCTAGAGGTTAATGAAGAAGGCTCAGAA
GCTGCTGTCTCAGGAATGATTGCAATTAGTAGGATGGCTGTGCTGTATCCTCAAGTTATTGTCGACC
ATCCATTTTCTTCTTATCAGAAACAGGAGAACTGGTACAATTCTATTCATGGGACGAGTCATGCATCC
TGAACAATGAACACAAGTGGACATGATTCGAAGAACTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC225640 representing NM_001122752
Red=Cloning site Green=Tags(s)


MAFLGLFSLLVLQSMATGATFPPEEAIADLSVNMYNRLRATGEDENILFSPLSIALAMGMELGAQGSTQE
 EIRHSMGYDSLKNGEEF SFLKEFSNMVTAKESQYVMKIANSLFVQNGFHVNEEF LQMMKKYFNAAVNHVD
 FSQNVAVANYINKWVENNTNLLVKDLVSPRDFDAATYLALINAVYFKGNWKSQFRPENTRTFSFTKDDDES
 EVQIPMMYQQGEFYGFSDGSNEAGGIYQVLEIPYEGDEISMMLVLSRQEVPLATLEPLVKAQLVEEWA
 NSYKKQKVEVYLPRTVEQEIDLKDVLLKALGITEIFIKDANLTGLYDNKEIFLSKAIHKSFLLEVNEEGSE
 AAASVGMIAISRMAVLYPQIVDHPFFFLIRNRRTGTILFMGRVMHPETMNTSGHDFEEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001122752

ORF Size:	1230 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001122752.1 , NP_001116224.1
RefSeq ORF:	1233 bp
Locus ID:	5274
UniProt ID:	Q99574
Cytogenetics:	3q26.1
Protein Families:	Druggable Genome, Secreted Protein
MW:	46.5 kDa
Gene Summary:	This gene encodes a member of the serpin superfamily of serine proteinase inhibitors. The protein is primarily secreted by axons in the brain, and preferentially reacts with and inhibits tissue-type plasminogen activator. It is thought to play a role in the regulation of axonal growth and the development of synaptic plasticity. Mutations in this gene result in familial encephalopathy with neuroserpin inclusion bodies (FENIB), which is a dominantly inherited form of familial encephalopathy and epilepsy characterized by the accumulation of mutant neuroserpin polymers. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]