

Product datasheet for RC225667L3

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc.

https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BAAT (NM_001127610) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: BAAT (NM_001127610) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: BAAT

Synonyms: BACAT; BACD1; BAT; HCHO

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Sgfl-Mlul

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC225667).

Sequence:

Restriction Sites:

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001127610

ORF Size: 1254 bp





BAAT (NM_001127610) Human Tagged Lenti ORF Clone - RC225667L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001127610.1</u>, <u>NP 001121082.1</u>

RefSeq Size: 3377 bp
RefSeq ORF: 1257 bp
Locus ID: 570

 UniProt ID:
 Q14032

 Cytogenetics:
 9q31.1

Protein Pathways: Biosynthesis of unsaturated fatty acids, Metabolic pathways, Primary bile acid biosynthesis,

Taurine and hypotaurine metabolism

MW: 46.3 kDa

Gene Summary: The protein encoded by this gene is a liver enzyme that catalyzes the transfer of C24 bile

acids from the acyl-CoA thioester to either glycine or taurine, the second step in the

formation of bile acid-amino acid conjugates. The bile acid conjugates then act as a detergent in the gastrointestinal tract, which enhances lipid and fat-soluble vitamin absorption. Defects in this gene are a cause of familial hypercholanemia (FHCA). Two transcript variants encoding

the same protein have been found for this gene. [provided by RefSeq, Jul 2008]