

Product datasheet for **RC225682**

MSF (SEPT9) (NM_001113492) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MSF (SEPT9) (NM_001113492) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SEPTIN9
Synonyms:	AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>NM_001113492 ORF sequence, RC225682 may differ due to SNPs. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCCCCCTGCCTCCAAGGTCCTCCGAGGTGCCACTGCCCTGCCACCGACGCGACCCCCAAGAGG
GTGGAGATCCAGATGCCAAGCCTGCTGAGGCGCCACCAGCCCCAGCCAGCCAGACCTTGGAGAAT
TCAGAGCCTGCCCTGTGTCTCAGCTGCAGAGCAGGCTGGAGCCCAAGCCCCAGCCCCCTGTGGCTGAG
GCTACACCCCGGAGCCAGGAGGCCACTGAGGCGGCTCCAGCTGCGTTGGCGACATGGCCGACACCCCC
AGAGATGCCGGGCTCAAGCAGGCGCCTGCATCACGGAACGAGAAGGCCCGGTGGACTTCGGCTACGTG
GGGATTGACTCCATCCTGGAGCAGATGCCCGGAAGGCCATGAAGCAGGGCTTCGAGTTCAACATCATG
GTGGTCGGGCAGAGCGGCTTGGGTAATCCACCTTAATCAACACCCTCTTCAAATCCAAAATCAGCCGG
AAGTCGGTGCAGCCACCTCAGAGGAGCGCATCCCCAAGACCATCGAGATCAAGTCCATCACGCACGAT
ATTGAGGAGAAAGGCGTCCGGATGAAGCTGACAGTGATTGACACACCAGGGTTCGGGGACCACATCAAC
AACGAGAACTGCTGGCAGCCCATCATGAAGTTCATCAATGACCAGTACGAGAAAATACCTGCAGGAGGAG
GTCAACATCAACCGCAAGAAGCGCATCCCGGACACCCGCTCCACTGCTGCCTCTACTTCATCCCGCC
ACCGGCCACTCCCTCAGGCCCTGGACATCGAGTTTATGAAACGCCTGAGCAAGTGGTCAACATCGTC
CCTGTCTCGCCAAGGCGGACACACTCACCTGGAGGAGAGGGTCCACTTCAAACAGCGGATCACCAGCA
GACCTGCTGCCAACGGCATCGACGTGTACCCCCAGAAGGAATTTGATGAGGACTCGGAGGACCGGCTG
GTGAACGAGAAGTCCGGGAGATGATCCCATTTGCTGTGGTGGCAGTGACCACGAGTACCAGGTCAAC
GGCAAGAGGATCCTTGGGAGGAAGACCAAGTGGGGTACCATCGAAGTTGAAAACACCACACTGTGAG
TTTGCCTACCTGCGGGACCTTCTCATCAGGACGCACATGCAGAACATCAAGGACATCACCAGCAGCATC
CACTTCGAGGCGTACCGTGTGAAGCGCCTCAACGAGGGCAGCAGCCCATGGCCAACGGCATGGAGGAG
AAGGAGCCAGAAGCCCCGGAGATG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RC225682
 Blue=ORF Red=Cloning site Green=Tag(s)


MEPPASKVPEVPTAPATDAAPKRVEIQMPKPAEAPTAPSPAQTLNENEPAPVSQLQSRLEPKPQPPVAE
 ATPRSQEATEAAPSCVGMADTPRDAGLKQAPASRNEKAPVDFGYVGIDSILEQMRKAMKQGFENIM
 VVGQSGLGKSTLINTLFKSKISRKSVQPTSEERIPKTIEIKSITHDIEEKGVRMKLTVIDTPGFGDHN
 NENCWQPIMKFINDQYEKYLQEEVNINRKKRIPDTRVHCLLYFIPATGHSLRPLDIEFMKRLSKVNVNIV
 PVIKADTLTLEERVHFKQRITADLLSNGIDVYPQKEFDEDEDLNEKFRMIPFAVVGSDHEVQVN
 GKRI LGRKTKWGTIEVENTTHCEFAYLRDLLIRTHMQNIKDITSSIHFEAYRVKRLNEGSSAMANGMEE
 KEPEAPEM
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001113492

ORF Size:	1266 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113492.2
RefSeq Size:	3953 bp
RefSeq ORF:	1269 bp
Locus ID:	10801
UniProt ID:	Q9UHD8
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome
MW:	47.5 kDa
Gene Summary:	This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2009]