

Product datasheet for **RC225703L3V**

NARGIL (NAA16) (NM_001110798) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NARGIL
Synonyms:	NARGIL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110798
ORF Size:	1287 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225703).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001110798.1](#)

RefSeq ORF: 1290 bp

Locus ID: 79612

UniProt ID: [Q6N069](#)

Cytogenetics: 13q14.11

Protein Families: Druggable Genome, Transcription Factors



MW: 50.4 kDa

Gene Summary: Auxillary subunit of the N-terminal acetyltransferase A (NatA) complex which displays alpha (N-terminal) acetyltransferase activity.[UniProtKB/Swiss-Prot Function]