

Product datasheet for **RC225704L3V**

ERLEC1 (NM_001127398) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ERLEC1 (NM_001127398) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ERLEC1
Synonyms:	C2orf30; CIM; CL24936; CL25084; HEL117; XTP3-B; XTP3TPB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001127398
ORF Size:	1287 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225704).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001127398.2 , NP_001120870.1
RefSeq Size:	2443 bp
RefSeq ORF:	1290 bp
Locus ID:	27248
UniProt ID:	Q96DZ1
Cytogenetics:	2p16.2
Protein Families:	Secreted Protein
MW:	48.9 kDa


[View online »](#)

Gene Summary:

This gene encodes a resident endoplasmic reticulum (ER) protein that functions in N-glycan recognition. This protein is thought to be involved in ER-associated degradation via its interaction with the membrane-associated ubiquitin ligase complex. It also functions as a regulator of multiple cellular stress-response pathways in a manner that promotes metastatic cell survival. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 21. [provided by RefSeq, Aug 2011]