

Product datasheet for **RC225706L4V**

GCNT7 (NM_080615) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GCNT7 (NM_080615) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GCNT7
Synonyms:	C20orf105; dj1153D9.2; gcnt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080615
ORF Size:	1290 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225706).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080615.1 , NP_542182.1
RefSeq ORF:	1293 bp
Locus ID:	140687
UniProt ID:	Q6ZNI0
Cytogenetics:	20q13.31
MW:	49.1 kDa
Gene Summary:	Glycosyltransferase.[UniProtKB/Swiss-Prot Function]


[View online »](#)