

Protein Sequence:	>RC225708 representing NM_001128925 Red=Cloning site Green=Tags(s) MDPAVALVLCLSLFLLSLWRQSSGRGRLPSGPTPLPIIGNILQLDVKDMSKSLTNFSKVYGPVFTVYFG LKPIIVLHGYEAVKEALIDHGEEFSGRGSFPVAEKVNKGLGILFSNGKRWKEIRRFCLMTLRNFGMGKRS IEDRVQEEARCLVEELRKTNASPCDPTFILGCAPCNVICSVIFHDRFDYKDQRFNLMEKFENLRILSS PWIQEKHNQQSEFTVESLIATVTDMFGAGTETTSTTLRYGLLLLLLKYPEVTAKVQEEIECVVGRNRSPCM QDRSHMPYTDVVHEIQRYIDLLPTNLPHAVTCDVKFKNYLIPKGTITIITSLTSVLHNDKEFPNPEMFDP GHFLDKSGNFKKSDYFMPFSAGKRMCMGEGLARMELFLFTTILQNFNLKSQVDPKIDITPIANAFGRV PPLYQLCFIPV TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001128925
ORF Size:	1293 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001128925.1
RefSeq ORF:	1296 bp
Locus ID:	1562
UniProt ID:	P33260
Cytogenetics:	10q23.33
MW:	48.6 kDa