

Product datasheet for **RC225719L4V**

PIP3E (IPCEF1) (NM_001130700) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PIP3E (IPCEF1) (NM_001130700) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PIP3E
Synonyms:	PIP3-E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001130700
ORF Size:	1314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225719).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130700.1 , NP_001124172.1
RefSeq Size:	6818 bp
RefSeq ORF:	1317 bp
Locus ID:	26034
UniProt ID:	Q8WWN9
Cytogenetics:	6q25.2
MW:	49.1 kDa



[View online »](#)

Gene Summary:

Enhances the promotion of guanine-nucleotide exchange by PSCD2 on ARF6 in a concentration-dependent manner.[UniProtKB/Swiss-Prot Function]