

Product datasheet for **RC225750**

ALDH3A1 (NM_001135167) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALDH3A1 (NM_001135167) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ALDH3A1
Synonyms:	ALDH3; ALDHIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC225750 representing NM_001135167
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCAAGATCAGCGAGGCCGTGAAGCGCGCCCGCGCCCTTCAGCTCGGGCAGGACCCGTCCGCTGC
 AGTTCCGGATCCAGCAGCTGGAGGCGCTGCAGCGCCTGATCCAGGAGCAGGAGCAGGAGCTGGTGGGCGC
 GCTGGCCGCAGACCTGCAACAAGATGAATGGAACGCTACTATGAGGAGGTGGTGTACGTCTAGAGGAG
 ATCGAGTACATGATCCAGAAGCTCCCTGAGTGGGCGCGGATGAGCCCGTGGAGAAGACGCCCCAGACTC
 AGCAGGACGAGCTCTACATCCACTCGGAGCCACTGGGCGTGGTCTCGTATTGGCACCTGGAACCTACC
 CTTCAACCTCACCATCCAGCCATGGTGGGCGCCATCGCTGCAGGGAACCTCAGTGGTCTCAAGCCCTCG
 GAGCTGAGTGAGAACATGGCGAGCTGCTGGCTACCATCATCCCCAGTACCTGGACAAGGATCTGTACC
 CAGTAATCAATGGGGGTGTCCCTGAGACCAGGAGCTGCTCAAGGAGAGGTTCCACCATATCCTGTACAC
 GGGCAGCACGGGGGTGGGAAGATCATCATGACGGCTGCTGCCAAGCACCTGACCCTGTACGCTGGAG
 CTGGGAGGGAAGAGTCCCTGCTACGTGGACAAGAACTGTGACCTGGACGTGGCCTGCCGACGCATCGCT
 GGGGAAATTCATGAACAGTGGCCAGACCTGCGTGGCCCTGACTACATCCTCTGTGACCCCTCGATCCA
 GAACAAATTTGTGGAGAAGCTCAAGAAGTCACTGAAAGAGTTCTACGGGGAAGATGCTAAGAAATCCCGG
 GACTATGGAAGAATCATTAGTCCCGGCACTTCCAGAGGGTGTGGGCTGATTGAGGGCCAGAAGGTGG
 CTTATGGGGGCACCGGGATGCCGCCACTCGCTACATAGCCCCACCATCCTCACGGACGTGGACCCCA
 GTCCCCGTGATGCAAGAGGAGATCTTCGGGCTGTGCTGCCATCGTGTGCGTGCAGCCTGGAGGAG
 GCCATCCAGTTCATCAACCAGCGTGAGAAGCCCTGGCCCTCTACATGTTCTCCAGCAACGACAAGGTGA
 TTAAGAAGATGATTGCAGAGACATCCAGTGGTGGGGTGGCGGCCAACGATGTCATCGTCCACATCACCT
 GCACTCTCTGCCCTTCGGGGCGTGGGAACAGCGGCATGGGATCCTACCATGGCAAGAAGAGCTTCGAG
 ACTTTCTCTACCGCCGCTCTTGCTGGTGGGCTCTGATGAATGATGAAGGCCTGAAGGTCAGATACC
 CCCCAGCCCGCCAAGATGACCCAGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC225750 representing NM_001135167
 Red=Cloning site Green=Tags(s)

MSKISEAVKRRAAFSSGRTRPLQFRIQQLEALQRLIQEQEQELVGALAADLHKNEWNAYEEVVYVLEE
 IEYMIQKLPEWAADPEVEKTPQTQDEL YIHSEPLGVVLVIGTWNYPFNLTIQPMVGAIAGNSVVLKPS
 ELSENMASLLATIIPQYLDKDLYPVINGGVPETTELLKERFDHILYTGSTGVGKIIMTAAAKHLTPVTL
 LGGKSPCYVDKNCDLVACRRIAWGKFMNSGQTCVAPDYILCDPSIQNQIVEKLLKSLKEFYGEDAKKSR
 DYGRII SARHFQRYMGLIEGQKVAYGGTDAATRYIAPTILTDVDPQSPVMQEEIFGPVLPIVCVRSLEE
 AIQFINQREKPLALYMFSSNDKVIKKMIAETSSGGVAANDVIVHITLHSLPFGGVSNGMGSYHGKKSFE
 TFSHRSCLVRPLMNDEGLKVRYPPSPAKMTQH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6150_f04.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001135167

ORF Size: 1359 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135167.1](#), [NP_001128639.1](#)

RefSeq ORF: 1362 bp

Locus ID: 218

UniProt ID: [P30838](#)

Cytogenetics: 17p11.2


Protein Families: Druggable Genome

Protein Pathways: Drug metabolism - cytochrome P450, Glycolysis / Gluconeogenesis, Histidine metabolism, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Phenylalanine metabolism, Tyrosine metabolism


MW: 50.2 kDa

Gene Summary: Aldehyde dehydrogenases oxidize various aldehydes to the corresponding acids. They are involved in the detoxification of alcohol-derived acetaldehyde and in the metabolism of corticosteroids, biogenic amines, neurotransmitters, and lipid peroxidation. The enzyme encoded by this gene forms a cytoplasmic homodimer that preferentially oxidizes aromatic and medium-chain (6 carbons or more) saturated and unsaturated aldehyde substrates. It is thought to promote resistance to UV and 4-hydroxy-2-nonenal-induced oxidative damage in the cornea. The gene is located within the Smith-Magenis syndrome region on chromosome 17. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2008]


Product images:



Circular map for RC225750



Western blot validation of overexpression lysate (Cat# [LY427570]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC225751] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified ALDH3A1 protein (Cat# [TP325750]). The protein was produced from HEK293T cells transfected with ALDH3A1 cDNA clone (Cat# RC225750) using MegaTran 2.0 (Cat# [TT210002]).