

Product datasheet for **RC225761**

RUNDC3B (NM_001134405) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	RUNDC3B
Synonyms:	RPIB9; RPIP9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC225761 representing NM_001134405
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCCTCCCGAGCCTGGGGGCTGAGCGGGATCCGCGCGGTGGCGCGGAGGCGGCAAGAAAAGCC
 TGAGCGCCCGCAATGCTGCGGTGGAGAGGAGGAACCTGATCACCCTGTGCAGGTTTTCTGTGAAGACCCT
 GATTGATCGGTCTTGCTTTGAGACAATTGATGATTCTTCTCCTGAATTTAACAATTTGCAGCTATTTTG
 GAACAGATTTTAAGCCACCGGCTAAAAGGTCAAGTAACCTGGTTTGGTTATGAAAGTCCTCGTAGCTTCT
 GGGACTATATCAGAGTGGCTTGCCGAAAGTTTACAGAATTGTATCTGCAGCATTGAAAAATATGAAAA
 TGTCAGTTCTTCTAGAGCTAAGGGTAGAGCCTGGATCAGAGTAGCACTCATGGAAAAACATTTATCTGAA
 TACATCTCTACAGCTCTGAGAGACTTCAAAACAACCAGGAGATTTTATGAAGATGGAGCAATTGTCTTG
 GTGAAGAAGCAATATGCTTGTGGCATGCTTCTAGGACTCAATGCTATTGATTTAGTTTCTGCCTAAA
 GGGAGAGGGGCTGGATGGCAGTTTTCTGCTGTAATAGACTATACCCATATTTGAAGTATATCCAAAGT
 TCTGATAGTATCAGCAGTGATGAGGAGGAGCTAAGGACTTTGGGAAGCAGTGGTAGCGAAAGCAGTACTC
 CAGAGAATGTCGGACCTCCTTTCTCATGGATGAGAACAGTTGGTTCAACAAGTGTAAGAGAGTTAAACA
 AAAGTATCAGCTTACCCTGGAACAGAAGGGTTACCTTGAAGAACTTACGACTTCGAGAGAACCACTA
 TCTGAATCTGTCTCCAGAATAAACTACTTCAAAGGATTGAAGATTCGATCTGGCTCATAAAGTGG
 AGAAGGAACAATTAGAATATATAATTGTGGAGCTTCAAGATCAGCTGACTGTGCTAAAGAATAATGATTT
 AAGATCGAGACAAGAGTTAACTGCCATCTCACCAACCAGTGGCCTTCTCCAGGAGCTCTGGATGTCAAT
 GCTGTTGCCTTGGATACGTTGCTTTACCGAAAACACAATAAACAGTGGTATGAGAAAAGTTATCAAAGTC
 TTGACCAGTTATCAGCAGAAGTAGCCTTTCTCAGACTTCACTAGATCCAGGCCAGTCACAAGAAGGAGA
 TGGAAACAAGACACATTAATGTAATGAGTGAAGGTAAAGGAAGATACTCCCTCATTACTTGGCCTCTGT
 GGATCTCTAACGTCAGTGGCAAGTTACAAGTCTCTAACAGCTTAAATCTAATGACTACCTTGCAAGTC
 CTACAACAGAGATGACAAGTCCAGGCTAACTCCATCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225761 representing NM_001134405
 Red=Cloning site Green=Tags(s)

MASRLGGLSGIRGGGGGGKKSLSARNAAVERRNLITVCRFSVKTLIDRSCFETIDDSSPEFNFAAIL
 EQILSHRLKGQVTWFGYESPRSFWDYIRVACRKVSQNCICSIENMENVSSRAKGRAWIRVALMEKHLSE
 YISTALRDFKTRRFYEDGAIVLGEEANMLAGMLLGLNAIDFSFCLKGEGLDGSFPAVIDYTPYLKYIQS
 SDSISSDEEELRTLGSSESSTPENVPFPLMDENSWFNKCKRVKQYQLTLEQKGYLEELLRLRENQL
 SESVSQNKILLQRIEDSLAHKLEKEQLEYIIVELQDQLTVLKNNDLRSRQELTAHLTNQWPSPGALDVN
 AVALDTLLYRKHNKQWYEKSYQSLDQLSAEVLSQTSLDPGQSQEGDGKQDTLNMSEGKEDTPSLLGLC
 GSLSVASYSKSLTSLKSNDYLASPTTEMTSPGLTPS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8061_b10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001134405

ORF Size: 1368 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001134405.2](#)


RefSeq ORF: 1371 bp

Locus ID: 154661
 UniProt ID: [Q96NLO](#)
 Cytogenetics: 7q21.12
 MW: 50.6 kDa

Product images:



Circular map for RC225761



Western blot validation of overexpression lysate (Cat# [LY427421]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225761 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).