

Product datasheet for **RC225773A1V**

Human **BIGM103 (SLC39A8) (NM_001135146) AAV Particle**

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	BIGM103
Synonyms:	BIGM103; CDG2N; LZT-Hs6; PP3105; ZIP8
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC225773 representing NM_001135146
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCCCCGGGTCGCGGGTGGCCGGGCTCCTGTTGCTGGCGGCCCGGCCCTCGGAGGAGTGCCGGAGG
GGCCAGGGCTAGCCTTCAGCGAGGATGTGCTGAGCGTGTTCGGCGCAATCTGAGCCTGTCGGCGGCCGA
GCTCCAGCACTTGCTGGAGCAGATGGGAGCCGCCTCCCGCTGGGCGTCCCGGAGCCTGGCCAGCTGCAC
TTCAACCAGTGTTAACTGCTGAAGAGATCTTTCCCTTCATGGCTTTTCAAATGCTACCCAAATAACCA
GCTCCAAATCTCTGTCTATCTGTCCAGCAGTCTTACAGCAATTGAACCTTCACCCATGTGAGGATCGGCC
CAAGCACAACCAAGACCAAGTCATTACAGAAGTTGGGGATATGGATTCCTGTCACTGACGATTATTAAT
CTGGCATCTCTCCTCGGATTGATTTTGACTCCACTGATAAAGAAATCTTATTTCCCAAAGATTTTGACCT
TTTTTGTTGGGGCTGGCTATTGGGACTCTTTTTTCAAATGCAATTTTCAAACCTATTCCAGAGGCATTTGG
ATTTGATCCCAAAGTCGACAGTTATGTTGAGAAGGCAGTTGCTGTGTTGGTGGATTTTACCTACTTTTC
TTTTTTGAAAGAATGCTAAAGATGTTATTAAGACATATGGTCAGAATGGTCATACCCACTTTGGAAATG
ATAACTTTGGTCTCAAGAAAAAACTCATCAACCTAAAGCATTACCTGCCATCAATGGTGTGACATGCTA
TGCAAATCCTGCTGTACAGAAGCTAATGGACATATCCATTTTGATAATGTCACTGTGGTATCTCTACAG
GATGGAAAAAAGAGCCAAGTTCATGTACCTGTTTGAAGGGGCCCAAACTGTCAGAAATAGGGACGATTG
CCTGGATGATAACGCTCTGCGATGCCCTCCACAATTTTCATCGATGGCCTGGCGATTGGGGCTTCCTGCAC
CTTGCTCTCCTTCAGGGACTCAGTACTTCCATAGCAATCCTATGTGAGGAGTTTCCCAACGAGTTAGGA
GACTTTGTGATCCTACTCAATGCAGGGATGAGCACTCGACAAGCCTTGCTATTCAACTTCCTTTCTGCAT
GTTCTGCTATGTTGGCTAGCTTTTGGCATTTTGGTGGGCAACAATTCGCTCCAATATTATATTTCG
ACTTGCTGGAGGCATGTTCTCTATATTTCTGTCGAGATATGTTTCCAGAGATGAATGATATGCTGAGA
GAAAAGGTAACGTGAAGAAAAACCGATTTACCTTCTTCATGATTGAGAATGCTGGAATGTTAACTGGAT
TCACAGCCATTCTACTCATTACCTTGATGTCAGGAGAAATCGAATTGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC225773 representing NM_001135146 Red=Cloning site Green=Tags(s) MAPGRAVAGLLLLAAAGLGGVAEGPGLAFSEDLVSVFGANLSLSAAQLQHLLQMGAAASRVGVPEPGQLH FNQCLTAEEIFSLHGFSNATQITSSKFSVICPAVLQQLNFHPCEDRPHKHKTRPSHSEVWGYGFLSVTIIN LASLLGLILTPLIKKSYFPKILTFVGLAIGTLFSNAIFQLIPEAFGDPKVDSEYVEKAVAVFGGFYLLF FFERMLKMLLKTYGQNGHTFGNDNFGPQEKTHQPKALPAINGVTCYANPAVTEANGHIHFDNVSVVSLQ DGKKEPSSCTCLKGPKLSEIGTIAWMITLCDALHNFIDGLAIGASCTLSLLQGLSTSIAILCEEFPHELG DFVILLNAGMSTRQALLFNFLSACSCYVGLAFGILVGNNFAPNIIIFALAGGMFLYISLADMFPEMNDMLR EKVTGRKTDFTFFMIQNAGMLTGFTAILLITLYAGEIELE TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001135146
ORF Size:	1380 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001135146.1
RefSeq ORF:	1383 bp
Locus ID:	64116
UniProt ID:	Q9C0K1
Cytogenetics:	4q24
MW:	49.5 kDa