

Product datasheet for **RC225774L3V**

FBRS (NM_001105079) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	FBRS
Synonyms:	FBS; FBS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001105079
ORF Size:	1380 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225774).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001105079.1](#), [NP_001098549.1](#)

RefSeq ORF: 2943 bp

Locus ID: 64319

UniProt ID: [Q9HAH7](#)

Cytogenetics: 16p11.2

Protein Families: Secreted Protein



MW: 48.2 kDa

Gene Summary: Fibrosin is a lymphokine secreted by activated lymphocytes that induces fibroblast proliferation (Prakash and Robbins, 1998 [PubMed 9809749]).[supplied by OMIM, Mar 2008]