

Product datasheet for **RC225790L3V**

GABRR3 (NM_001105580) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GABRR3 (NM_001105580) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GABRR3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001105580
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225790).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001105580.1
RefSeq ORF:	1404 bp
Locus ID:	200959
UniProt ID:	A8MPY1
Cytogenetics:	3q11.2
Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors
MW:	54.1 kDa



[View online »](#)

Gene Summary:

The neurotransmitter gamma-aminobutyric acid (GABA) functions in the central nervous system to regulate synaptic transmission of neurons. This gene encodes one of three related subunits, which combine as homo- or hetero-pentamers to form GABA(C) receptors. In humans, some individuals contain a single-base polymorphism (dbSNP rs832032) that is predicted to inactivate the gene product. [provided by RefSeq, Jan 2012]