

## Product datasheet for **RC225803L3V**

### PRAMEF17 (NM\_001099851) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Symbol:</b>                   | PRAMEF17                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Mammalian Cell Selection:</b> | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Vector:</b>                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Tag:</b>                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>ACCN:</b>                     | NM_001099851                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ORF Size:</b>                 | 1422 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>ORF Nucleotide Sequence:</b>  | The ORF insert of this clone is exactly the same as(RC225803).                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>OTI Disclaimer:</b>           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| <b>OTI Annotation:</b>           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq:</b>                   | <a href="#">NM_001099851.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>               | 1425 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>                 | 391004                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>UniProt ID:</b>               | <a href="#">Q5VTA0</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Cytogenetics:</b>             | 1p36.21                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>MW:</b>                       | 54.6 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                           |

