

Product datasheet for **RC225809**

WDSUB1 (NM_001128213) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	WDSUB1
Synonyms:	UBOX6; WDSAM1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC225809 representing NM_001128213
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGTGAACTGATTACACATTAGCTGATCATGGTGACGATGTCAACTGCTGTGCCTTCTCCTTTCC
TCTTGCTACTTGCTCCTTGGACAAAACAATTCGCCTGTACTCGTTACGTGACTTTACTGAACTGCCACA
TTCTCCATTGAAGTTTCATACCTATGCTGTCCACTGCTGCTGTTTCTCCCTTCAGGACATATTTTGGCA
TCGTGTTCAACAGATGGTACCCTGTCTTATGGAATGAAATGGACAGATGCTGGCAGTGATGGAAC
AGCCTAGTGGCAGCCCTGTGAGGGTTTCCAGTTCCTCCAGACTCCACGTGTTTGGCATCAGGGGCAGC
TGATGGAACGTGGTGTGTTTGTGGAATGCACAGTCATACAAATTATATAGATGTGGTAGTGTAAAGATGGC
TCCTTGGCGCATGTGCATTTCTCTAATGGAAGCTTCTTGTCACTGGCTCCTCATGTGGTGATTAA
CAGTGTGGGATGATAAATGAGGTGTCTGCATAGTGAAGACACATGATCTTGAATTACCTGCTGCGA
TTTTTCTTACAGCCAGTTTCTGATGGAGAACAGGTCTTCAGTTTTTTCGACTGGCATCATGTGGTCAG
GATTGCCAAGTCAAAATTTGATTGTTTCTTTACCCATATCTTAGGTTTTGAATTAATAAAGTA
CACTGAGTGGGCACTGTGCTCCTGTTCTGGCTTGTGCTTTTCCCATGATGGGCAGATGCTAGTCTCAGG
GTCAGTGGATAAGTCTGTCAATAGTATATGATACTAATACTGAGAATATACTTCACACATTGACTCAGCAC
ACCAGGTATGTCAACTTGTGCTTTTGCACCTAATACCTTTTACTTGCTACTGGTTCAATGGACAAAA
CAGTGAACATCTGGCAATTTGACCTGGAAACACTTTGCCAAGCAAGGCGCACAGAACATCAGCTGAAGCA
ATTTACCGAAGATTGGTCAGAGGAGGATGTCTCAACATGGCTTTGTGCACAAGATTTAAAGATCTTGTT
GGTATTTTCAAGATGAATAACATTGATGGAAAAGAACTGTTGAATCTTACAAAAGAAAGTCTGGCTGATG
ATTTGAAAATGAATCTCTAGGACTGCGTAGTAAAGTCTGAGGAAAATTGAAGAGCTCAGGACCAAGGT
TAAATCCCTTTCTCAGGAATTCCTGATGAATTTATATGTCCAATAACTAGAGAACTTATGAAAGATCCG
GTCATCGCATCAGATGGCTATTCATATGAAAAGGAAGCAATGGAATTTGGATCAGCAAAAAGAAACGTA
CAAGTCCCATGACAAATCTGTTCTTCTCAGCGGTACTTACACCAATAGGACTCTGAAAATGGCCAT
CAATAGATGGCTGGAGACACACCAAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225809 representing NM_001128213
Red=Cloning site Green=Tags(s)


MVKLIHTLADHGDDVNCCAFSFSLLATCSLDKTIIRLYSLRDFTELPHSPLKFHTYAVHCCCFSPSGHILA
SCSTDGTTVLWNTENGQMLAVMEQPSGSPVRVCQFSPDSTCLASGAADGTVVWNAQSYKLYRCGSVKDG
SLAACAFSPNGSFFVTGSSCGDLTVWDDKMRLHSEKAHDLGITCCDFSSQPVSDGEQGLQFRLASCGQ
DCQVKIWIWVSFTHILGFELKYKSTLSGHCAPVLACAFSHDGMQLVSGSVDKSVIVYDNTNENILHTLTQH
TRYVTTCAFPNTLLLATGSMCKTVNIWQFDLETLCQARRTEHQLKQFTEWSEEDVSTWLCAQDLKDLV
GIFKMNNIDGKELLNLTKESLADDLKIESLGLRSKVLRIEELRTKVKSLSSGIPDEFICPITRELMKDP
VIASDGYSEKEAMENWISKKKRTSPMTNLVLPASVLTNPRTLKMAINRWLETHQK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6598_h12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001128213

ORF Size: 1428 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001128213.2](#)

RefSeq ORF: 1431 bp

Locus ID: 151525
Cytogenetics: 2q24.2
Protein Families: Druggable Genome
MW: 52.7 kDa

Product images:

Circular map for RC225809

Western validation with an anti-DDK antibody
* L: Control HEK293 lysate R: Over-expression lysate