

Product datasheet for **RC225856L3V**

WDR26 (NM_001115113) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	WDR26
Synonyms:	CDW2; GID7; MIP2; SKDEAS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001115113
ORF Size:	1494 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225856).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001115113.1](#), [NP_001108585.1](#)

RefSeq ORF: 1938 bp

Locus ID: 80232

UniProt ID: [Q9H7D7](#)

Cytogenetics: 1q42.11-q42.12

MW: 56.94 kDa



Gene Summary:

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Two transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]