

Product datasheet for **RC225861L3V**

PHACS (ACCS) (NM_001127219) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	PHACS
Synonyms:	ACS; PHACS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001127219
ORF Size:	1503 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225861).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001127219.1
RefSeq ORF:	1506 bp
Locus ID:	84680
UniProt ID:	Q96QU6
Cytogenetics:	11p11.2
MW:	57.1 kDa



Gene Summary:

Does not catalyze the synthesis of 1-aminocyclopropane-1-carboxylate but is capable of catalyzing the deamination of L-vinylglycine.[UniProtKB/Swiss-Prot Function]