

Product datasheet for **RC225867**

NELF (NSMF) (NM_001130971) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NELF (NSMF) (NM_001130971) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NELF
Synonyms:	HH9; NELF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC225867 representing NM_001130971
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCGCCGCCCTCCAGGAGGAGGGCGCTGAGGAGCGAGGCCATGTCCTCGGTGGCGGCCAAAGTGC
 GAGCAGCCCGAGCGTTTGGAGAGTACCTGTCCCAGAGTACCCTGAGAACCGCAACGGCGAGATCACCT
 GCTGGCTGATGCCTACTCTGGCCACGACGGTCCCGGAGATGCAGCCGGCCCCCAGAACAAGCGCCGC
 CTGTCCCTCGTCTCCAACGGCTGCTACGAGGGCAGCCTCTCAGAGGAGCCAGCATTAGGAAGCCCGCAG
 GCGAGGGCCCTCAGCCTCGAGTGTACACCATCTCTGGGGAGCCTGCCCTGCTGCCAGCCCTGAGGCGGA
 GGCCATTGAGCTGGCGGTGGTGAAGGGGCGGGCAGCGGCACCCTCACCATCACAGCCAGCCCTGCGC
 GCCAGCCCTGGTGGCAGCCGGGAGGACGTGAGCAGGCCCTGCCAGAGCTGGGCGGCAGCCGCCAGGGCT
 CCAAGGAGTGCCCGGATGTGCCAGCTGGCTCCTGGCCCCACCCTCGGGCCTTTGGGTGGACCAGCC
 ACCTCTGCCTGAGACCTCCGGTCCCGCAAGAAGCTGGAGAGGATGTACAGCGTTGACCGTGTGTCTGAC
 GACATCCCTATTGCTACCTGGTTCCCAAGGAAAATCTTTTCAGCTTCCAGACAGCAACCACAACATATGC
 AAGCGAATTCGCAACACCTGCGCATGGTCGGCAGCCGGAGGGTGAAGGCCAGACGTTTCGCTGAGCG
 GCGCGAGCGGAGCTTACGCCGTCTGGAGCGACCCACCCCATGAAAGCCGACACTTCCACGACTCC
 CGAGACAGCAGTGACCTGCAGAGCTCCACTGCACGCTGGACGAGGCCTTCGAGGACCTGGACTGGGACA
 CTGAGAAGGGCCTGGAGGCTGTGGCCTGCGACACCGAAGGCTTCGTGCCACCAAAGGTATGCTCATTTC
 CTCAAGGTGCCAAGGCTGAGTACATCCCACATATCCCGCCGGGATGACCCCTCCATCATCCCCATC
 CTCTACGACCATGAGCACGCAACCTTCGAGGACATCCTTGAGGAGATAGAGAGGAAGCTGAACGTTACC
 ACAAGGGAGCCAAGATCTGGAATAATGCTGATTTTCTGCCAGGGAGTCTGGACACCTCTATCTCCTCAA
 GAACAAGGTGGCCACCTTTGCCAAAGTGGAGAAGGAAGGAGACATGATTCACTTCTGGAAGCGGCTGAGC
 CGCCTGATGAGCAAAGTGAACCCAGAGCCGAACGTCATCCACATCATGGGCTGCTACATTCTGGGAACC
 CCAATGGAGAGAAGCTGTTCCAGAACCTCAGGACCCTCATGACTCCTTATAGGGTACCTTCGAGTCACC
 CCTGGAGCTCTCAGCCCAAGGAAGCAGATGATCGAGACGTACTTTGACTTCCGGTTGTATCGCCTGTGG
 AAGAGCCGCCAGCACTCGAAGCTGCTGGACTTTGACGACGTCCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC225867 representing NM_001130971
 Red=Cloning site Green=Tags(s)

MGAAASRRRALRSEAMSSVAAKVRAARAFGEYLSQSHPENRNGADHLLADAYSGHDGSPMQPAPQNKRR
 LSLVSNCGYEGSLSEEPSIRKPAGEGPQPRVYTIISGEPALLPSPEAEIELAVVKGRRQRHPHHHSQPLR
 ASPGGSREDVSRPCQSWAGSRQGSKECPGCAQLAPGPTPRAFGLDQPPLPETSRRRKKLERMYSVDRVSD
 DIPIRTWFPKENLFSFQTATTTMQANFRKHLRMVGSRRVKAQTFERRERSFSRSWSDPTPMKADTSHDS
 RDSSDLQSSHCTLDEAFEDLDWDEKLEAVACDTEGFVPPKVMLISSKVPKAEYIPTIIRRDDPSIPI
 LYDHEHATFEDILEEIERKLNIVYHKGAKIWKMLIFCQGGPGHLYLLKNKVATFAKVEKEEDMIHFWKRLS
 RLMSKVNPEPNVIHIMGCYILGNPNGEKLFQNLRTLMTPLYRVTFESPLELSAQGKQMIETYDFRLYRLW
 KSRQHSKLLDFDDVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001130971

ORF Size: 1515 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130971.1](#), [NP_001124443.1](#)

RefSeq ORF: 1518 bp

Locus ID: 26012

UniProt ID: [Q6X4W1](#)


Cytogenetics: 9q34.3

Protein Families: Druggable Genome

MW: 57 kDa

Gene Summary: The protein encoded by this gene is involved in guidance of olfactory axon projections and migration of luteinizing hormone-releasing hormone neurons. Defects in this gene are a cause of idiopathic hypogonadotropic hypogonadism (IHH). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2010]

Product images:



Circular map for RC225867