

Product datasheet for **RC225896**

LRRC48 (DRC3) (NM_001130090) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	LRRC48
Synonyms:	CFAPI34; LRRC48
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RC225896 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACCAGCCGTGCAACTCGATGGAGCCGAGGGTGATGGACGATGACATGCTCAAGCTGGCCGTCGGGG
 ACCAGGGCCCCCAGGAGGAGGCCGGGCAGCTGGCCAAGCAGGAGGGCATCCTCTTCAAGGATGTCCTGTC
 CCTGCAGCTGGACTTTCGGAACATCTCCGCATAGACAACCTCTGGCAGTTTGAGAAGCTGAGGAAGCTG
 CAGCTGGACAATAACATCATTGAGAAGTCGAGGGCCTGGAGAACCTCGCACACCTGGTCTGGCTGGATC
 TGTCTTTCAACAACATTGAGACCATCGAGGGCTGGACACACTGGTGAACCTGGAGGACCTGAGCTTGT
 CAACAACCGGATCTCCAAGATCGACTCCCTGGACGCCCTCGTCAAGCTGCAGGTGTTGTCGCTGGGCAAC
 AACCGGATTGACAACATGATGAACATCATCTACCTTCGGCGTTCAAGTGCCTGCGGACGCTCAGCCTCT
 CTAGGAACCTATCTCTGAGGCAGAGGATTACAAGATGTTTCATCTGTGCCTACCTTCTGACCTCATGTA
 CCTGGACTACTGGCGATTGATGACCACACAAAAAGCTTGCAGGAGCTAAGCACCAGTACAGCATCGAC
 GAGCTGAAGCACCAGGAGAACCTGATGCAGGCCAGCTGGAGGACGAGCAGGCGCAGCGGGAGGAGCTAG
 AGAAGCACAAGACTGCGTTTGTGGAACACCTGAATGGCTCCTTCTGTTTGACAGCATGTACGCTGAGGA
 CTCAGAGGGCAACAATCTGTCTACCTGCTGCTGGTGTGCGGTGAGCTCCTTGAGACCTACAAGGACAAGTT
 GTCATCATCTGCGTGAATATTTTGTGATGGCTGAAACAGCAGGAGAAGCGGAAAAACAGAGCTTGACA
 CCTTCAGTGAATGTGTCGCTGAGGCCATCCAGGAAAAACAGGAGCAGGGCAACGCAAGATTGCCAAATT
 CGAGGAGAAGCACTTGTGAGTTTAAAGTGCCATTGAGAGGAGTTGGAAGTGCCTAACATTGAGAAGATG
 ATCCTAGAATGCAGTGTGACATCAGTGAGTTGTTGATGCGCTCATGACGCTGGAGATGCAGCTGGTGG
 AGCAGCTGGAGGAGACTATAACATGTTTGAAGGAACATTGTTGACATGGTAGGACTGTTTATCGAAAA
 TGTCCAAAGCCTGATGGCTCAGTGCCGGGACCTGGAGAATCACCACCACGAGAAGCTCCTGGAGATCTCT
 ATCAGCACCTGGAGAAGATTGTGAGGGCGACCTGGACGAGGACCTGCCTAACGACCTGCGCGCGCTTT
 TTGTCGATAAAGATACGATTGTTAATGCTGTCGGGGCATCGCACGACATCCACCTCCTGAAGATTGACAA
 TCGAGAAGATGAGCTGGTGACCAGAATCAACTCTTGGTGACACGTTTAAAGACAGGATTACAAGGAT
 GAGATCATGAGGAACCGAAGCGCGTGAAGGAGATCAATCAGTACATCGACCACATGCAGAGCGAACTGG
 ACAACCTGGAATGTGGCGACATCCTAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225896 protein sequence
 Red=Cloning site Green=Tags(s)


MNQPCNSMEPRVMDMLKLA VGDQGPQEEAGQLAKQEGILFKDVL SLQLDFRNILRIDNLWQFENLRKL
 QLDNNIEKIEGLENL AHLVWDL SFNNIETIEGLDTLVNLEDL SLFNNRISKIDSLDALVKLQVLSLGN
 NRIDNMMNIYLRFKCLRTL SL SRNPISEAEDYKMFICAYLPDL MYLDYWRIDDHTKKLAEAKHQYSID
 ELKHQENLMQAQLEDEQAQREELEKHKTAFVEHLNGSFLFDSMYAEDSEGNLSYLPVGELLEITYKDKF
 VIICVNI FEYGLKQEKRTKELDTFSECVREAIQENQE QGKRKIAKFEEKHLSSL SAIREELELPNIEKM
 ILECSADISELFDALMTLEMLVEQLEETINMFERNIVDMVGLFIENVQSLMAQCRDLENHHHEKLLLEIS
 ISTLEKIVEGDLDEDL PNDLRALFVDKDTIVNAV GASHDIHLKIDNREDELVTRINSWCTRLIDRIHKD
 EIMRNRKRKVEINQYIDHMQSELDNLECGDILD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6349_all.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001130090

ORF Size: 1569 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22µm filter is required.


RefSeq: NM_001130090.1, NP_001123562.1

RefSeq Size: 2183 bp


RefSeq ORF: 1572 bp
Locus ID: 83450
UniProt ID: [Q9H069](#)
Cytogenetics: 17p11.2
MW: 61.1 kDa

Gene Summary: Component of the nexin-dynein regulatory complex (N-DRC) a key regulator of ciliary/flagellar motility which maintains the alignment and integrity of the distal axoneme and regulates microtubule sliding in motile axonemes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC225896



Western blot validation of overexpression lysate (Cat# [LY427158]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225896 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).